


Fourier 300

- Console
Service Manual
Version 003



Copyright © by Bruker Corporation

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form, or by any means without the prior consent of the publisher. Product names used are trademarks or registered trademarks of their respective holders.

This manual was written by

Diether Maier and Frank Decker

© October 28, 2014 Bruker Corporation

P/N: Z31980

DWG: Z4D11520

For further technical assistance for this product, please do not hesitate to contact your nearest BRUKER dealer or contact us directly at:

Bruker Corporation

Am Silberstreifen

D-76287 Rheinstetten

Germany

Phone: +49 721 5161 0

FAX: +49 721 5171 01

E-mail: nmr-support@bruker.de

Internet: www.bruker.com

Contents

1	About This Manual	7
1.1	Policy Statement	7
1.2	Symbols and Conventions	7
1.3	Font and Format Conventions	9
2	Introduction	11
2.1	Concept	11
2.2	Limitation of Liability	11
2.3	Before you Begin	12
2.4	Minimum Qualifications for Service Personnel	12
2.5	The Bruker Service	13
2.6	Transport to Manufacturer	13
3	Safety	15
3.1	General	15
3.2	Personnel Requirements	15
3.2.1	Qualifications	15
3.2.2	Unauthorized Persons	16
3.2.3	Instruction	16
3.3	Personal Protective Equipment	16
3.4	Basic Dangers	16
3.4.1	General Workplace Dangers	17
3.4.2	Dangers from Electric Power	18
3.4.3	Mechanical Dangers	19
3.4.4	Dangers from Gases Under Pressure	19
3.4.5	Dangers from Magnetic Fields	20
3.4.6	Dangers Due to High or Low Temperatures	21
3.4.7	Danger from Chemical Substances	21
3.5	Environmental Protection	22
3.6	Signage	23
3.7	Spare Parts	23
4	Tool Requirements	25
5	Design and Function	27
5.1	Overview	27
5.2	Brief Description	27
5.2.1	Internal Console Component Layout	29
5.3	Unit Descriptions	34
5.3.1	Intelligent Devices with Ethernet	34
5.3.1.1	Main Board	34
5.3.1.2	Lock System	35
5.3.1.3	Variable Temperature Unit (VTU)	36
5.3.2	Intelligent Devices with Can-Bus	36
5.3.2.1	Shim Control	36

5.3.2.2	Peripheral Control Board	37
5.3.2.3	GAB/2 Backplane	37
5.3.3	Intelligent Devices with LVDS	37
5.3.4	Devices without Interface	37
5.3.4.1	Transmitters	37
5.3.4.2	Preamplifiers	38
5.3.4.3	Lock Switch	38
5.3.4.4	Power Supplies	40
5.3.4.5	Ventilators	41
5.4	Connections	42
5.5	Accessories	42
6	Installation.....	43
6.1	Unpacking	43
6.2	Wiring Up	43
6.3	Adjusting Air Pressure	44
6.4	Switch on for the First Time	44
6.5	Checking that the Console has Booted	44
6.6	Performing cf on a Fourier	45
6.6.1	Setting the Spectrometer Frequency	46
6.6.2	ETH Addresses and tty Settings	46
6.6.3	Cortab and Power Check	47
6.6.4	Editing the Nuclei List	47
6.6.5	RF Routing	49
6.7	Additional Configurations	50
6.7.1	Installation of Pulse Programs and Parameter Sets	50
6.7.2	Editing Solvents	52
6.7.3	Configuring the Probe	52
6.7.4	Configuring the Lock Parameters	54
6.7.5	Configuring the Transmitter and Preamplifier Pre-delays	56
6.8	Setting the Magnet to the Right Field	56
6.9	The Lock System	57
6.9.1	Working with the Lock System	57
6.9.2	Finding a Lock	57
6.10	Converting the Parameter Set to the Actual BF1	63
6.11	Tuning and Matching the Probe	64
6.11.1	General Information about the Fourier Probe	64
6.12	The Gradshim Command	64
6.13	Installation Acceptance	66
6.14	Main Differences between Fourier & Avance Systems	68
7	Configuration.....	69
7.1	Network Settings	69
7.1.1	Resetting the CPU	70
7.2	Firmware Update	71
7.2.1	Ethernet Devices (Main Board, Lock Board, VTU)	71
7.2.2	CAN-Bus Devices	72
7.3	FPGA Update	74

7.3.1	Lock Control and VTU	74
7.3.2	Main Board	75
7.4	Device Specific Settings & Tools	76
7.4.1	Main Board	77
7.4.1.1	BIS EEPROM Settings.....	77
7.4.1.2	Board Temperature Control	79
7.4.1.3	On-Board Voltages	80
7.4.1.4	Pulse Programmer Memory	80
7.4.2	Shim Control	81
7.4.3	GAB/2 Interface	84
7.4.4	Lock Board	85
7.4.5	VTU	92
7.4.5.1	Changing PID Parameters	93
7.4.5.2	Setting the Heater Protection Threshold.....	94
7.4.5.3	Calibrating the Thermocouple.....	95
7.4.6	Peripheral Control	95
7.4.6.1	Setting the Lift Parameters	95
7.4.6.2	Setting the Spin Parameters	96
7.4.6.3	Setting the Spin Safety Threshold	97
7.4.6.4	Calibrating and Measuring the Helium Level	97
7.4.6.5	Monitoring the Console and Shim Stack Temperature	99
7.5	Additional Test Applications	100
8	Pulse Programming.....	105
8.1	General Principle	105
8.2	Structure of Pulse Programs	105
8.3	Syntax and Examples	106
9	Wiring Schedules	111
10	Maintenance and Repair	119
11	Cleaning and Disposal	121
11.1	Cleaning	121
11.2	Disposal	121
11.3	Safety	121
11.4	Dismantling	122
11.5	Disposal Instructions	122
12	Technical Specifications.....	123
13	Contact	127
	List of Figures.....	129
	List of Tables	131
	Index	133

1 About This Manual

This manual is intended to be a reference guide for operators and service technicians. It provides detailed information about the user level maintenance and service, as well as overall use of the Bruker device.

This manual enables safe and efficient handling of the device.

This manual is an integral part of the device, and must be kept in close proximity to the device where it is permanently accessible to personnel. In addition, instructions concerning labor protection laws, operator regulations tools and supplies must be available and adhered to.

This manual is planned for use by the service engineer in the installation of the device or accessory, and for the user in performing existing experiments and for the setup of new experiments.

This manual is an integral part of the device, and must be kept in close proximity to the device where it is permanently accessible to personnel. In addition, instructions concerning labor protection laws, operator regulations tools and supplies must be available and adhered to.

Before starting any work, personnel must read the manual thoroughly and understand its contents. Compliance with all specified safety and operating instructions, as well as local accident prevention regulations, are vital to ensure safe operation.

The figures shown in this manual are designed to be general and informative and may not represent the specific Bruker model, component or software/firmware version you are working with. Options and accessories may or may not be illustrated in each figure.

1.1 Policy Statement

It is the policy of Bruker to improve products as new techniques and components become available. Bruker reserves the right to change specifications at any time.

Every effort has been made to avoid errors in text and figure presentation in this publication. In order to produce useful and appropriate documentation, we welcome your comments on this publication. Support engineers are advised to regularly check with Bruker for updated information.

Bruker is committed to providing customers with inventive, high quality products and services that are environmentally sound.

1.2 Symbols and Conventions

Safety instructions in this manual are marked with symbols. The safety instructions are introduced using indicative words which express the extent of the hazard.

In order to avoid accidents, personal injury or damage to property, always observe safety instructions and proceed with care.

DANGER



DANGER indicates a hazardous situation which, if not avoided, will result in death or serious injury.

This is the consequence of not following the warning.

1. This is the safety condition.
 - ▶ This is the safety instruction.

WARNING



WARNING indicates a hazardous situation, which, if not avoided, could result in death or serious injury.

This is the consequence of not following the warning.

1. This is the safety condition.
 - ▶ This is the safety instruction.

CAUTION



CAUTION indicates a hazardous situation, which, if not avoided, may result in minor or moderate injury.

This is the consequence of not following the warning.

1. This is the safety condition.
 - ▶ This is the safety instruction.

NOTICE

NOTICE indicates a property damage message.

This is the consequence of not following the notice.

1. This is a safety condition.
 - ▶ This is a safety instruction.

SAFETY INSTRUCTIONS

SAFETY INSTRUCTIONS are used for control flow and shutdowns in the event of an error or emergency.

This is the consequence of not following the safety instructions.

1. This is a safety condition.
 - ▶ This is a safety instruction.



This symbol highlights useful tips and recommendations as well as information designed to ensure efficient and smooth operation.

1.3 Font and Format Conventions

Type of Information	Font	Examples
Shell Command, Commands, “All what you can enter”	Arial bold	Type or enter fromjdx zg
Button, Tab, Pane and Menu Names “All what you can click”	Arial bold, initial letters capitalized	Use the Export To File button. Click OK . Click Processing...
Windows, Dialog Windows, Pop-up Windows Names	Arial, initial letters capitalized	The Stacked Plot Edit dialog will be displayed.
Path, File, Dataset and Experiment Names Data Path Variables Bruker Trademarks Table Column Names Field Names (within Dialog Windows)	Arial Italics	<i>\$tshome/exp/stan/nmr/</i> <i>lists</i> <i>expno, procno,</i> <i>IconNMR™,</i> <i>TopSpin™,</i> <i>XWIN-NMR™</i>
Parameters	Arial in Capital Letters	VCLIST
Program Code Pulse and AU Program Names Macros Functions Arguments Variables	Courier	go=2 au_zgte edmac CalcExpTime() XAU(prog, arg) disk2, user2
AU Macro	Courier in Capital Letters	REX PNO

Table 1.1: Font and Format Conventions

2 Introduction

This service manual is intended to be used by trained Bruker service staff only. The manual is a service and installation guide for the console part of the Fourier 300 system. Please refer to the individual probe and magnet manuals for service related information on these devices.

2.1 Concept

The Fourier 300 is intended for use as an FT-NMR spectrometer for structure verification in basic research, routine NMR measurements in chemical laboratories, and for educational purposes in the field of NMR. It delivers powerful performance at extremely compact size. With its new robust probe technology and a unique push-button, power on/off concept, ease of siting and handling is guaranteed. Researchers have access to numerous pre-defined 1D and 2D experiments and interactive, automated processing tools help to transfer spectroscopic data into a corresponding report.

The Fourier 300 console consisting of:

- Fourier 300 spectrometer.
- EASY NMR probe (dual probe).
- Actively shielded superconducting magnet.
- BOSS1 shim system.



Figure 2.1: The Fourier 300 System

2.2 Limitation of Liability

All specifications and instructions in this manual have been compiled taking account of applicable standards and regulations, the current state of technology and the experience and insights we have gained over the years.

The manufacturer accepts no liability for damage due to:

- Failure to observe this manual
- Improper use
- Deployment of untrained personnel

- Unauthorized modifications
- Technical modifications
- Use of unauthorized spare parts

The actual scope of supply may differ from the explanations and depictions in this manual in the case of special designs, take-up of additional ordering options, or as a result of the latest technical modifications.

The undertakings agreed in the supply contract as well as the manufacturer's Terms and Conditions and Terms of Delivery and the legal regulations applicable at the time of conclusion of the contract shall apply.

2.3 Before you Begin

This service manual contains information and safety information that are necessary for the installation and servicing of the device.

All maintenance and repairs are to be accomplished using the information in this manual. At the same time references over general maintenance and care from the User Manual are also to be followed.

Consider all safety references!

Be sure that no NMR measurement is in progress during the installation as it might be disturbed.

Information for ordering spare parts is available in the spare parts section from the Bruker Service Center (see [Contact](#) [127]).

2.4 Minimum Qualifications for Service Personnel

Type of Task	Personnel	Training and Experience
Transportation	No special requirements.	No special.
Installation	Bruker certified personnel only.	Technically skilled, with a good knowledge of the application field.
Routine Use Daily Maintenance	Appropriately certified and experienced personnel, familiar with use of computers and automation in general	Laboratory technicians or equivalent. Training is usually done in-house. Familiar with MS Windows® environment.
Setup and optimization of program	Bruker certified personnel only.	Experienced laboratory technician. High degree of knowledge of the relevant application field.
Preventive Maintenance	Bruker certified personnel only.	Technically skilled with a basic understanding of the application.
Servicing	Bruker certified personnel only.	Background and experience in electronics/mechanics with computer knowledge.

Table 2.1: Overview Installation and Operation Requirements for Personnel

2.5 The Bruker Service

Our customer service division is available to provide technical information. See [Contact \[▶ 127\]](#) for details.

In addition, our employees are always interested in acquiring new information and experience gained from practical application; such information and experience may help improve our products.

2.6 Transport to Manufacturer

When the device must be returned to the manufacturer for a major repair, use the original packaging for transportation.

Include a good description of the problem.

3 Safety

This section provides an overview of all the main safety aspects involved in ensuring optimal personnel protection and safe and smooth operation.

Non-compliance with the action guidelines and safety instructions contained in this manual may result in serious hazards.

3.1 General

Before you start any repair inside of the device, be aware of the high 230/115V voltages. Even if these voltages are protected by security features to avoid any physical contact, it is still possible that the voltage sources can be unintentionally touched with a tool, object, etc.

Therefore, always check if you really need the power supply to be switched on during your work. Otherwise turn the device off and disconnect the power cable from the wall socket to the device. Safeguard that no one is able to re-power the system without your approval.

3.2 Personnel Requirements



Only trained Bruker personnel are allowed to mount, retrofit, repair, adjust and dismantle the unit!

3.2.1 Qualifications

This manual specifies the personnel qualifications required for the different areas of work, listed below:


Laboratory Personnel


Laboratory personnel are health care professionals, technicians, and assistants staffing a research or health care facility where specimens are grown, tested, or evaluated and the results of such measures are recorded. Laboratory personnel are able to carry out assigned work and to recognize and prevent possible dangers self-reliant due to their professional training, knowledge and experience as well as profound knowledge of applicable regulations.

The workforce must only consist of persons who can be expected to carry out their work reliably. Persons with impaired reactions due to, for example, the consumption of drugs, alcohol, or medication are prohibited from carrying out work on the device.

When selecting personnel, the age-related and occupation-related regulations governing the usage location must be observed.

3.2.2 Unauthorized Persons

 **WARNING**



Risk to life for unauthorized personnel due to hazards in the danger and working zone!

Unauthorized personnel who do not meet the requirements described in this manual will not be familiar with the dangers in the working zone. Therefore, unauthorized persons face the risk of serious injury or death.

- ▶ Unauthorized persons must be kept away from the danger and working zone.
- ▶ If in doubt, address the persons in question and ask them to leave the danger and working zone.
- ▶ Cease work while unauthorized persons are in the danger and working zone.

3.2.3 Instruction

Personnel must receive regular instruction from the owner. The instruction must be documented to facilitate improved verification.

Date	Name	Type of Instruction	Instruction Provided By	Signature

3.3 Personal Protective Equipment

Personal protective equipment is used to protect the personnel from dangers which could affect their safety or health while working.

Personnel must wear personal protective equipment while carrying out the different operations at and with the device.

This equipment will be defined by the head of the laboratory. Always comply with the instructions governing personal protective equipment posted in the work area.

3.4 Basic Dangers

The following section specifies residual risks which may result from using the device and have been established by means of a risk assessment.

In order to minimize health hazards and avoid dangerous situations, follow the safety instructions specified here as well as in the following chapters of this manual.

3.4.1 General Workplace Dangers

WARNING



Danger to life from nonfunctional safety devices!

If safety devices are not functioning or are disabled, there is a danger of serious injury or death.

- ▶ Check that all safety devices are fully functional and correctly installed before starting work.
- ▶ Never disable or bypass safety devices.
- ▶ Ensure that all safety devices are always accessible.

CAUTION



Danger of injury from tripping over dirt and scattered objects!

Dirt and scattered objects may cause people to slip or trip, resulting in personal injuries.

- ▶ Always keep the work area clean.
- ▶ Remove objects which are no longer required from the work area and particularly from the floor.
- ▶ Indicate unavoidable hazards using marking tape.

CAUTION



Accident hazard from falling from ladder!

It is possible to fall from a ladder when it is used to reach the device on some magnets.

- ▶ Do not use a ladder.
- ▶ Use an approved platform to reach the device on the magnet.
- ▶ Wear non-slip shoes.

NOTICE

Material damage due to a software error!

Samples or the device may be damaged due to a software error causing malfunction of the control system. Users may also be shocked by abrupt malfunction or unexpected system start.

- ▶ Dummy samples must be used during installation and service.
- ▶ Personnel should be alerted to unexpected malfunctions.

NOTICE

Material damage hazard due to impacting the magnet!

Impacting the magnet may result in a quench.

- ▶ Mount the device carefully on the magnet.
- ▶ Avoid banging the magnet during installation and operation, e.g. when replacing the sample carousel.

NOTICE

Material damage due to the use of genuine samples during installation and maintenance!

Using genuine samples during installation and maintenance may result in material damage.

- ▶ Use only dummy samples during installation and maintenance.

3.4.2 Dangers from Electric Power

DANGER

Danger to life from stored charges!

Electric charges may be stored in electrical components even after the system has been switched off and disconnected from the power supply. Contact with these components may result in serious or fatal injury.

- ▶ Before working on the specified components, ensure that they have been completely disconnected from the power supply.
- ▶ Allow 10 minutes to elapse in order to ensure that the internal capacitors have been fully discharged.



WARNING

Danger of injury from electrical shock!

A life threatening shock may result when the housing is open during operation.

- ▶ Only qualified personnel should open the housing.
- ▶ Disconnect the device from the electrical power supply before opening the device. Use a voltmeter to verify that the device is not under power!
- ▶ Be sure that the power supply cannot be reconnected without notice.



WARNING

Danger to life from residual electrostatic potentials!

Friction between material being conveyed may result in significant development of electrostatic potential. Therefore, contact with parts immediately following the conveying operation may be life-threatening.

- ▶ Potential equalisation must be ensured before making contact with parts, unless such equalisation is provided by the customer.



Electrostatic discharge from friction may occur, resulting in an electric spark and loud bang. Use ESD flooring and wear ESD shoes.

⚠ WARNING**Danger to life from contact voltage!**

Absent or faulty protective earth conductor may result in contact voltage. This may pose a risk of injury or death.

- ▶ Before the initial commissioning of the device, connect the main power supply to the socket and verify the complete functionality of the protective earth conductor.

3.4.3 Mechanical Dangers**⚠ CAUTION****Accident hazard from movement of mechanical parts!**

The fingers or hand may be pinched due to movement of mechanical parts.

- ▶ Shut off the device before accessing.

⚠ CAUTION**Accident and material damage hazard from falling objects!**

Equipment may fall down during assembly, retrofitting, or dismantling. This may result in personal injury or equipment damage.

- ▶ If necessary, assemble/disassemble the device in multiple parts.
- ▶ Use a platform with railings instead of a ladder to reach the assembly area.
- ▶ Avoid working over the head. When this can not be avoided, wear a protective hard hat.
- ▶ Follow the mounting instructions in the installation manual.

3.4.4 Dangers from Gases Under Pressure**⚠ WARNING****Danger of injury due to movements caused by stored pneumatic forces!**

Pneumatically driven components may move unexpectedly due to stored residual forces, causing serious injuries.

- ▶ Work on the pneumatics system must only be carried out by trained pneumatics technicians.
- ▶ Before starting work on the pneumatics system, ensure that it has been completely depressurized. The pressure accumulator must be completely relieved.

WARNING



Accident hazard from asphyxiation!

A break in the pneumatic hose may result in the uncontrolled exit of nitrogen into the laboratory.

- ▶ An oxygen warning device should be present in the laboratory if the device is operated with nitrogen.
- ▶ Note that leakage from the main supply line cannot be stopped by the device!

Overpressure Valve

The high pressure system includes an overpressure valve which safely reduces the excess pressure in the event of inadmissible pressure conditions developing as a result of faulty operation, component failure or other irregular events.

3.4.5 Dangers from Magnetic Fields

WARNING



Danger to life from strong magnetic fields!

Strong magnetic fields may cause serious injuries or death and significant damage to property.

- ▶ Persons fitted with heart pacemakers must be kept away from the device. The functionality of the heart pacemaker could be compromised.
- ▶ Persons with metal implants must be kept away from the device. Implants may heat up or be subject to magnetic attraction.
- ▶ Ferromagnetic materials, tools, and electromagnets must be kept away from the magnetic source. Such materials could be subject to magnetic attraction and may fly around the room, injuring or killing people. Minimum distance 3 meters.
- ▶ Remove magnetic items (jewelry, watches, pens etc.) before carrying out maintenance work.
- ▶ Keep electronic equipment away from the magnetic source. Such equipment could be damaged.
- ▶ Keep storage media, credit cards etc. away from the magnetic source. Data could be erased.



The magnetic field of the device does not cause any personal injuries or property damage. For further information see the manual of the magnet used.

3.4.6 Dangers Due to High or Low Temperatures



CAUTION

Accident hazard from contact with hot or cold surfaces on the unit!

Contact with the hot or cold surfaces of the unit may result in serious burns.

- ▶ Do not touch device parts of cooled or heated units.
- ▶ Do not use damaged samples.
- ▶ After removing a sample or cassette allow it to cool or thaw before coming in contact.

NOTICE

Material damage hazard from overflow of cryogenics.

Material damage may result from the overflow of cryogenics.

- ▶ Turn off the device during magnet servicing.
- ▶ Cover the device with a protective cover to avoid contact with cold gases.
- ▶ Be sure to use sufficient transfer line and Teflon evacuation hose for nitrogen and helium refills based on recommendations in the magnet manual.
- ▶ After refilling cryogenics some parts of the magnet may be icy. Be sure to remove the ice to avoid its melting onto the device.

3.4.7 Danger from Chemical Substances



DANGER

Danger of injury from glassware breakage!

Broken glassware may cause minor injuries or material damage, but may also result in a life threatening situation if hazardous substances are used.

- ▶ If glassware breaks, refer to the corresponding precautions and cleaning/disinfection instructions.
- ▶ Wear protective equipment.
- ▶ Perform all tasks with the glassware carefully.
- ▶ Before carrying out any maintenance work, remove the samples and use dummy samples if necessary.
- ▶ Strictly observe the correct sample adjustment, i.e. the maximum sample height.
- ▶ Always transport the glassware with the cover, if applicable. Never turn the glassware upside down or on its side.
- ▶ **The laboratory supervisor is responsible for:**
 - ⇒ Establishing and enforcing standard sample handling and cleaning procedures.
 - ⇒ Establishing and enforcing the use of protective clothing and equipment.
 - ⇒ Training laboratory personnel.
 - ⇒ Preparing an emergency plan.



WARNING

Danger of injury from vapor formation!

During the work process, vapors may form which cause serious injury if inhaled.

- ▶ Only install the device in a well-ventilated room or ensure that an extractor is fitted.

NOTICE

Material damage hazard from material contact with NMR solvents!

Material damage may result when the device comes in contact with NMR solvents.

- ▶ Follow instructions provided in the manual for correct handling of solvents.
- ▶ Follow the sensor cleaning procedures described in this manual.
- ▶ If surface damage should occur, contact Bruker for repair of damaged parts.

NOTICE

Material damage hazard from heavy samples!

Samples may be damaged due to incorrect sample lift pressure adjustment.

- ▶ Adjustment is valid only for 1 sample configuration and weight.
- ▶ Personnel must be trained.

3.5 Environmental Protection

NOTICE

Danger to the environment from incorrect handling of pollutants!

Incorrect handling of pollutants, particularly incorrect waste disposal, may cause serious damage to the environment.

- ▶ Always observe local environmental regulations regarding handling and disposal of pollutants.
- ▶ Take the appropriate actions immediately if pollutants escape accidentally into the environment. If in doubt, inform the responsible municipal authorities about the damage and ask about the appropriate actions to be taken.

The following pollutants are used:

Nitrogen gas	Nitrogen gas may cause suffocation at high concentrations. Disposal of the empty gas cylinders must be performed by a specialist disposal company.
Helium inert gas	Helium inert gas may cause suffocation at high concentrations. Disposal of the empty gas cylinders must be performed by a specialist disposal company.

- Coolants** When released, coolants develop decomposition products which are hazardous to the environment. Maximum care and caution are required when handling coolants. Always observe the safety data sheet issued by the manufacturer. Ensure that personnel handling coolants are regularly informed about potential dangers and are instructed in the safe handling of coolants.
- Cleaning liquids** Cleaning liquids incorporating solvents contain toxic sub-stances. They must not be allowed to escape into the environment. Disposal must be carried out by a specialist disposal company.

3.6 Signage

The following symbols and information signs can be found in the work area. They refer to their immediate surroundings.



The identification and placement of warning labels are included in the manual. The laboratory supervisor is responsible for ensuring that all the warning labels are maintained in their proper place any time that the device is used.

Electrical Voltage



Only qualified electricians are permitted to work in a work room marked by this sign. Unauthorized persons must not enter the workplaces thus marked and must not open the marked cabinet.

Danger Spot



Warning indicating a danger spot in work rooms.
The warning label may be ordered using Bruker Part Number 67470.

3.7 Spare Parts

Loss of Guarantee

The use of non-approved spare parts will invalidate the manufacturer's guarantee.
Purchase spare parts from authorized dealers or directly from the manufacturer. See Contact for manufacturer's address.

4 Tool Requirements

All tools required for the installation are provided in the accessory case.



Be sure to use non-magnetic tools, especially when working on unshielded magnets.

5 Design and Function

5.1 Overview



Figure 5.1: The Fourier 300 Console

5.2 Brief Description

+The Fourier 300 Console consists of the following main components:

- [Intelligent Devices with Ethernet \[34\]](#)
- [Intelligent Devices with Can-Bus \[36\]](#)
- [Intelligent Devices with LVDS \[37\]](#)
- [Devices without Interface \[37\]](#)

The following block diagram shows how all the components work together:




Figure 5.2: Fourier 300 Console Overview

5.2.1 Internal Console Component Layout




Figure 5.3: Console: Left Side

1	ON/OFF Switch	7	1H Transmitter
2	Shim Control Boards (2)	8	28 Vdc PSU
3	Lock Control	9	Lock Switch
4	Lock Preamp	10	Peripheral Control
5	Lock Atten.	11	C13 Preamp
6	C13 Transmitter	12	2 x 24 Vdc Fans




Figure 5.4: Console: Right Side

1	+5V, +/- 12V, 24V PSU	4	Variable Temp Unit (VTU)
2	1H Preamp	5	Main Board
3	Pneumatic Valves	6	Main Board CPU




Figure 5.5: Console: Front View

1	Ethernet Hub	4	GAB/2 Backplane
2	Air Inlet (air filter not shown)	5	2H Lock Preamp
3	GAB/2	6	Power ON/OFF Control Board



Figure 5.6: Location of the Console Air Filter Tray

1	Air Filter Tray Remove the tray in order to clean the air filter.
---	----------------------------------------------------------------------




Figure 5.7: Console Input and Output Connections

1	To/From 2H Probe	10	To Sample Spin
2	Mains IN	11	To Sample Lift
3	To/From C13 Probe	12	N2 Gas Flow to Probe

4	To BOSS1 Shim System	13	To Heater
5	From Helium Sensor	14	From Thermocouple
6	From Sample Detector	15	Ethernet Out (e.g. to Sample Changer)
7	To/From 1H Probe	16	Ethernet To/From PC
8	To Gradient Coil	17	Accessory Mains Out

5.3 Unit Descriptions

5.3.1 Intelligent Devices with Ethernet

The intelligent devices with Ethernet connectors are connected to the Front End PC via the built-in hub. These devices include:

- Main Board (Pulse Programmer, 1H Modulator, C13-Modulator, Receiver, Digital Filter)
- Lock System (Digital 2H Lock, including Pre-amplifier and power attenuator)
- VTU (Variable Temperature Unit for N2 Sample Temperature Control)

5.3.1.1 Main Board

The Main Board contains the Pulse generation for two RF channels, the Receiver Channel including digital filtering and Gradient Channel with LVDS Interface.

The Pulse Programmer is preloaded with a complete 1D- or 2D-Experiment. Once started, the entire experiment is run from the FPGA without interaction from the CPU. There is an external trigger that can be used to trigger individual scans.

The receiver uses a direct sampling technique to digitize and down convert the RF signal. The FPGA contains a DSP core which performs the digital filtering down the final required spectral bandwidth. The data from each single FID is stored in the FPGAs internal memory. Accumulation of FIDs is performed by the on board Arm CPU, which then transfers the data to the front end whenever it is requested.

The main board has a fixed transmit but a variable receiver routing.

DDS2 is equivalent to FCU2 in the signal routing table and is dedicated to 1H.




Figure 5.8: DDS2

DDS1 is equivalent to FCU1 in the signal routing table and is dedicated to C13.




Figure 5.9: DDS1

The Main Board is also the Can-Bus Master for the shim Control and the GAB/2 Backplane.

The main board receives only 28Vdc from the central 28V PSU and generates all the other voltages it requires on board.

5.3.1.2 Lock System

The Lock System comprises three unit that are mounted together on one module:

- Lock-Board (transmit and receive)
- Lock Power Attenuator
- 2H Preamplifier

The lock board generates continuous small flip angle pulses at nominal 1 kHz rate. The lock transmit power can be varied in the range from 0dBm to -30 dBm using the separate digital attenuator. The receiver used a 50 MHz ADC to digitally down-convert the FID using the same principle as on the main board. The thus obtained real part of the FID is then used for the lock monitor, whereas the imaginary part is used to track the field, using a digital P+I control.

The lock operates on a field locking principle, using the H0 Coil of the BOSS1 Shim System. The total range of the lock is around 20 kHz (as measured on 1H). This is controlled by a 12-bit coarse DAC for range selection and a 12-bit fine DAC which does the actual field regulation. The resolution of the fine DAC is about 0.025 Hz per DAC Step. In case of magnet drift, the system has an **auto-range** function. i.e. when the fine DAC gets to below 500 or above 3500, the coarse DAC is incremented or decremented by one step, to bring the fine DAC back nearer the center of its operating region.

The two DAC voltages are fed to the shim board where they are scaled and added. A current driver provides the driver current for the H0 coil.

The Lock board receives only 28Vdc from the central 28V PSU and generates all the other voltages it requires on board.

5.3.1.3 Variable Temperature Unit (VTU)

The VTU is used to control the sample temperature. The heater and thermocouple connections are fully compatible with the BVT3000. It supports Thermocouple Types T and K but unlike the BVT3000 it has an internal electronic cold junction compensation, which means that ambient temperature variations do not affect the reading. The gas flow regulation is done using 4 digital valves which allow 16 gas flow settings from 0 to 2000/h. The recommended gas flow for the Fourier “one for all” probe is between 670 and 800 l/h. For setting PID values please refer to the VTU web server description.

5.3.2 Intelligent Devices with Can-Bus

- 2 x Shim Boards (12 Channels each + H0 Control) – Slave of Main Board
- Peripheral Control Board (Spin, Lift, Helium Level Control, Sample Status) – Slave of VTU
- GAB/2 Backplane (Interface to GAB/2 Board and -24V power Supply) Slave of Main Board

5.3.2.1 Shim Control

The Shim Control comprises two shim boards, each of which is capable of controlling up to 12 channels and a H0 coil. The two boards are cascaded and controlled from the main board CPU via CAN-Bus. Each channel has a 16-bit DAC and can deliver up to 300 mA. The shim board has an ADC with a multiplexer which can measure the current flowing through each coil for diagnostic purposes. Please refer to the web server description for details.

The shim currents are fed to the peripheral control board via 4 ribbon cables (2-per board) where they are converted to the standard BOSS1 Connector.

The 20 shim terms are divided up between two shim boards:




Figure 5.10: Two Shim Boards

The shim control boards are identical in build. The definition as to whether a board is shim board #1 or #2 is done via a jumper connector which is plugged into ST11. If this connector is present, it becomes board #2. This is also shown by LED LD7 which is then lit.

Since shim board #1 also controls the H0 current, it has a connection to the Lock Board, where the lock control voltage is generated.

The shim boards receive only 28Vdc from the central 28V PSU and generate all the other voltages they requires on board.

5.3.2.2 Peripheral Control Board

The peripheral control board takes care of spin, lift, sample control, helium level measurement and also routes the individual shim cables to a standards BOSS1-Connector. It is controlled via CAN-BUS from the **VTU**. It also contains temperature sensors to monitors the console and shim stack temperatures. Spin and Lift are controlled using proportional pneumatic valves that require a PWM (pulse width modulation) to set the required airflow. Please refer to the VTU web server description on how to set up the optimal air flow settings for spin and lift, as well as how to calibrate the helium level sensor.

5.3.2.3 GAB/2 Backplane

The GAB Backplane provides the GAB/2 Gradient Amplifier with all the necessary infrastructure to initialize and calibrate the GAB/2 gradient amplifier, i.e. it performs the function normally carries out by the ELCB Board. It also contains a DC/DC converter to provide the -24VDC which are otherwise not present in the console. It is controlled via CAN-Bus from the main board and has an additional slot for an optional ELCB Board for debug purposes.

The actual control of the gradient information is done via the LVDS interface from the main board.

5.3.3 Intelligent Devices with LVDS

- GAB/2 Gradient Amplifier

5.3.4 Devices without Interface

- 2 x Transmitters (1H & C13-Channel)
- 1H Preamplifier
- C13 Preamplifier
- Lock Switch (to route C13 Transmitter and Preamp to 2H Lock Coil for gradient shimming purposes)
- $\pm 5V$, 12V & 24V Power Supply (used for VTU, Preamps and Ethernet switch)
- 28 V Power Supply (used for everything else)
- 2 x Ventilators

5.3.4.1 Transmitters

There are two transmitters which are identical in build but have a fixed routing so that the upper one is dedicated to 1H (providing about 50W at 300 MHz) and the other is dedicated to C13 (providing about 150W at 75 MHz). The transmitters receiver their RF and gating inputs from the main board and have an RF-Output monitor which can be routed to the main board receiver for diagnostic purposes.

The Transmitters receives only 28Vdc from the central 28V PSU. The gating inputs on both transmitters are active high TTL.

5.3.4.2 Preamplifiers

There are three dedicated preamplifiers. The 13C (75 MHz) and 2H (46 MHz) Preamp are of the same type. The 1H (300 MHz) preamp is unique. All preamps have a gain of ca. 40dB.

Each amp has an RF in, RF out and two logic inputs. "Gate" is used to actively gate the preamp for low power pulses that don't have enough voltage to gate themselves. This is a coaxial input and is driven by one of the service pulses on the main board.

Tune/Blank is a dual purpose input. It switches the amplifier to *tune* mode, which measures the reflective power and is used for the wobble function. Since this mode also bypasses the 40 dB gain stage, it also blanks the preamp. The software always blanks the preamp that is not being used for observing, to avoid reflected power from decoupling pulses on the other channel getting to the receiver. The tune/gate signal is generated by the main board and is fed to the lock switch via a ribbon cable. From there it routes the tune/gate signal to the individual preamps (see wiring schedule in appendix). LEDs indicate the status of tune/blank mode (see below).

The preamps do not have a separate protection input as they are self-protecting. They receive +12V as a single supply voltage.

Function	Sub-D9	1H	13C / 2H
Tune/Blank	6	Hi	Hi
Gate	SMB	Low	Hi
GND	3;7;8	0V	0V
Vsupply	4;5;9	12V	12V

Table 5.1: Preamplifiers

5.3.4.3 Lock Switch

The Lock Switch is required in order to perform gradient shimming on Deuterium for samples that do not contain water. It reroutes the transmit and receive signals from the C13 channel to the 2H coil. For diagnostic purposes 2H observe spectroscopy can be carried out by selecting the 2H observe nucleus.




Figure 5.11: Lock Switch

The lock switch signal is generated by the main board and is fed to the lock switch via a ribbon cable. An audible click can be heard when it is activated. A LED also indicates the status of the lock switch - see the figure below.




Figure 5.12: Lock Switch Block Diagram




Figure 5.13: Lock Switch connections & LEDES

5.3.4.4 Power Supplies

The Fourier300 Console has two AC power supplies:

PSU 1: Input 90 – 250 Vac, Output: +5V, ±12V and +24Vdc (+standby voltage)

- VTU
- Fans
- Ethernet Hub
- Preamplifiers

PSU 2: Input 90 – 250 Vac, Output: +28Vdc, 20 A

- Main Board
- Lock Board
- Shim Boards
- Transmitters
- GAB/Backplane

When the mains switch on the rear panel is switched on, the console first goes into stand by which is indicated by a blinking ON/OFF switch. Pressing the ON-OFF Switch then provides both PSU1 and PSU2 via inhibit signals provided by the Power ON/OFF Control Board which is located close to the ON/OFF switch.

PSU1 is fused with 4A slow blowing fuses. The power consumption of the console is approx. 300 W.

For details of how the power is distributed to the units see the wiring schedule in the appendix.

5.3.4.5 Ventilators

Two 24Vdc Ventilators at the rear panel suck out warm air from the console. The ventilators have serial resistors to bring their actual operating voltage down to about 18V, where they operate more silently, prolong ventilator and filter lifetime and still create sufficient air flow.

The fresh air intake is in the base plate near the front of the console. The air filter can be removed for cleaning services without opening the console.



Figure 5.14: Ventilator

5.4 Connections




Figure 5.15: Connections

5.5 Accessories

The Fourier has the following accessories / options:

- Sample Express Lite Sample Changer
- BSCU N2 Cooling Unit
- Magnet vibration Isolation

6 Installation

6.1 Unpacking



WARNING

Do not lift the Fourier up and out of the carton!

Lifting the Fourier up and out of the carton may cause injury and damage the device.

- ▶ Simply slide the carton sleeve up and over the console, and then slide the console off the pallet.

NOTICE

Material damage hazard due to incorrect unpacking procedure!

Lifting the Fourier console out of the carton may damage the device.

- ▶ Simply slide the carton sleeve up and over the console, and then slide the console off the pallet.
- ▶ The Console has wheels and can be moved around freely using the handle at the back.
- ▶ Retain packaging in case the console has to be sent back for repair/exchange later.

6.2 Wiring Up

Please refer to the diagram in chapter [Wiring Schedules](#) [▶ 111].

External Filters

Some older models have three external filters (Z13383, Z14331, Z00116).

- Z13383 is connected to 1H Connector the rear of the console.
- Z14331 is connected to C13 Connector the rear of the console.
- Z00116 is connected to C13 Connector of the probe.

Newer models just have two external filters (E1451420, E1451430)

- E1451420 is connected to 1H Connector the rear of the console.
- E1451430 is connected to C13 Connector the rear of the console.

6.3 Adjusting Air Pressure

Input pressure to the Fourier 300 console:

- Minimum: 0.4 MPa (4 bar/58 psi)
- Maximum: 0.7 MPa (7 bar/100 psi)

Adjust the air pressure so that the manometer in the console reads between 3.5 and 4 bar. That way the default spin and lift parameters will get you started.

6.4 Switch on for the First Time

The console has three modes:

- OFF (Blue light OFF), switch at rear panel set to OFF.
- STAND BY (Blue light BLINKING), switch at rear panel set to OFF, ON/OFF Button at top of console not yet pressed.
- ON (Blue light ON).



It is IMPORTANT to switch the console on BEFORE switching on the PC, since the console contains the Ethernet hub, and the DHCP Server will not start if it doesn't find a network.

6.5 Checking that the Console has Booted

This procedure takes 30-60 seconds, provided the DHCP Server is running. There is no visual indication on the console, but the easiest way to check is to see if the three main Ethernet devices (Main Board, Lock and VTU) can be reached via ping or a web browser. Their addresses will usually be in the range:

149.236.99.25x

If you are not sure the DHCP Server is running, you can check by opening

Control Panel -> Administrative Tools -> Services -> Bruker DHCP Server




Figure 6.1: Checking the DHCP Server

If you are not sure what addresses have been assigned, you can check by opening the following file

Bruker\Diskless\dchpd.leases

```

lease 149.236.99.254 {
  starts 4 2010/08/19 11:08:55;
  ends 2 2038/01/19 03:14:06;
  tstp 2 2038/01/19 03:14:06;
  binding state active;
  next binding state free;
  hardware ethernet 00:00:ad:11:4f:12;
  uid "\001\000\000\255\0210\022";
}
lease 149.236.99.253 {
  starts 4 2010/08/19 11:08:58;
  ends 2 2038/01/19 03:14:06;
  tstp 2 2038/01/19 03:14:06;
  binding state active;
  next binding state free;
  hardware ethernet 00:00:ad:13:34:12;
  uid "\001\000\000\255\0234\022";
}
lease 149.236.99.252 {
  starts 4 2010/08/19 11:09:01;
  ends 2 2038/01/19 03:14:06;
  tstp 2 2038/01/19 03:14:06;
  binding state active;
  next binding state free;
  hardware ethernet 00:00:ad:12:85:12;
  uid "\001\000\000\255\022\205\022";
}

```

Once you are sure that all three devices have booted, you can start TopSpin. Normally your PC will come with **cf** having already been carried out. If you get an error message, that one or more devices cannot be found, you may have to repeat **cf**.

6.6 Performing cf on a Fourier

TopSpin 3.1 is the first software to be used with the Avance system and the Fourier system. All changes to TopSpin regarding the new Fourier console are implemented in the same software tree. That means that a Fourier spectrometer can be controlled from any computer running TopSpin 3.1.

As with the Avance systems, the Fourier must be configured in TopSpin using the **cf** command. The NMR administrator password is required to run **cf**.

The first screen allows to edit an already existing configuration or to create a new one.



Figure 6.2: cf

6.6.1 Setting the Spectrometer Frequency

The next screen asks for the configuration name and spectrometer frequency (1H) in MHz. It is usual for Avance300 systems to set it to 300.13MHz. The field of the magnet is close to the corresponding value and the correct field value is achieved by adding a B0 shim current to it. The Fourier console works the same way. But in order to keep the B0 shim current at a low value, the frequency should be set to a value close to the one which originates from the actual field of the respective magnet.

At this point, there is a minor inconsistency which can't be solved. The value of the magnetic field is checked using the NMR console to detect a proton or Deuterium resonance line. That means that it is necessary to run **cf** twice. The first run can be done with the default 1H frequency of 300.13MHz. Without using the lock, one has to search for the 1H spectrum of a suitable sample (Sucrose or Doped Water, use **zg** with a pulse length of about 5us and a spectrum width of >200kHz). The frequency at which the spectrum is detected, is close enough to the correct value. Choose a rounded value next to it and re-run **cf** with this value.

In the example shown here, the magnet's field resulted in a proton frequency of roughly 300.178MHz. This could be rounded up to 300.18 or down to 300.17, e.g., to have numbers being easy to memorize.

6.6.2 ETH Addresses and tty Settings

During the **cf** run, some information is shown about the detected hardware. In addition to the Fourier main board, only the lock system, the VTU and the Lock switch are part of the basic instrument. The Ethernet addresses of the Lock and the VTU are set by the DHCP server and don't need to be set manually. A missing entry shows that the respective board couldn't be detected. This may indicate a problem of the board or the DHCP server.

The Lock switch has to be set to amplifier 1, which is the X amplifier.

The preamplifiers and power amplifiers are build-in and are directly controlled by the Fourier main board.

Currently, none of the Avance accessories is used with the Fourier system. This may be done in the future, though.

6.6.3 Cortab and Power Check

Cortab files for the console are already measured as part of the production. So, please enable power check for probe protection.







Figure 6.3: POWCHK

6.6.4 Editing the Nuclei List

It might be a good idea to edit the nuclei list. The Fourier is a dedicated 1H, 13C and 2H instrument. So, only these nuclei are needed. This is not mandatory but a lot of selection windows are much easier to scroll if the nuclei list has been cut short.




When **cf** has completed its device search, it comes up with a page of information showing all the devices it found.




6.6.5 RF Routing

The routing window is rather easy. The green lines indicating a fixed hardware wiring have to be drawn as shown in the next screenshot.




The receiver routing is done automatically and is not shown by default. It can be made visible by clicking on **show receiver routing** and must look like this.



The configuration ends with a summary of all the hardware and settings.

6.7 Additional Configurations

After finishing, **cf** automatically opens a window from which additional configurations can be started.



6.7.1 Installation of Pulse Programs and Parameter Sets

Clicking on **Expinstall** starts the installation and conversion of pulse programs and datasets. Choose **Installation for Spectrometer**.

The next window will show an entry **Fourier** after finalizing TopSpin 3.1. click on it to install Fourier specific pulse programs and parameter sets.








Figure 6.4: expinstall

On the following window, one can choose which files should be installed: pulse programs, AU programs, CPD programs, shape files, etc. followed by the selection window for printer and plotter.

The last window lets one select the digitizer and default acquisition mode.

Follow the next few screen shots using default settings:





6.7.2 Editing Solvents


The list showing the common solvents in solution NMR is displayed and can be edited in the same way as the nucleus list.

6.7.3 Configuring the Probe

The command **edhead** is used to make the probe specific parameters known to the spectrometer. Modern NMR probes use a special memory called PICS to store these data. This PICS memory can be read by the HPPR preamplifier system via a cable connection. This is not available on the Fourier system. Therefore, all parameters must be either be set manually (see below) or you should copy the probe.ph file, which is delivered with the probe to the directory



conf/instr/FOURIER300


before running EDHEAD. Then select the probe you have installed. The probe which comes with the Fourier system is a DUL 13C/2H-1H Z-Grad probe.




Click on **Edit Probe Parameters** to define the probe. All values necessary for that are delivered on a spec sheet which can be found in the probe case.

The most important ones are the shown here:






The last step of probe definition is the signal routing to the preamplifiers which looks like the figure below:



6.7.4 Configuring the Lock Parameters


edlock gives access to the lock parameters being used for different solvents. The lock parameters for the Fourier are different from those for the Avance BSMS. Whereas the solvent related parameters as DISTANCE (in ppm) or FIELD-CORRECTION (in Hz) are unchanged, the hardware related parameters as LOCK POWER and the PID values are different. The Fourier comes with its own list which is installed by **cf**.



The screenshot shows the Edlock software window with the 'Solvents' tab selected. The window title is 'Edlock' and the menu bar includes 'Solvents', 'Edit', 'BSMS', 'Import', and 'Help'. Below the menu bar are tabs for 'Solvents', 'Lock', 'Signal', and 'Properties'. The main area contains a table with the following data:

△ Solvent	Lock Power	Loop Gain	Loop Time	Loop Filter	Lock Phase
Acetone	-29	0.4	0.4	4	-1
C6D6	-20	0.4	0.4	4	-1
C6D6+Dioxane	-20	0.4	0.4	4	-1
CD2Cl2	-20	0.4	0.4	4	-1
CDCl3	-12	0.4	0.4	1000	-1
D2O	-20	0.4	0.4	4	-1
Dioxane	-20	0.4	0.4	4	-1
DMSO	-20	0.4	0.4	4	-1
EtOD	-20	0.4	0.4	4	-1
H2O+D2O	-12	0.4	0.4	1000	-1
MeOD	-12	0.4	0.4	4	-1
Tol	-28	0.4	0.4	1000	-1


At the bottom right of the window are two buttons: 'Restore default' and 'Close'.



The screenshot shows the 'Edit solvent' dialog box for 'Acetone'. The title bar reads 'Edit solvent' and the window title is 'Edit solvent "Acetone"'. The dialog contains the following fields:


- Probe name: Fourier
- Probe description: Default probe
- Lock phase: -1
- Lock power: -29
- Loop filter: 4
- Loop gain: 0.4
- Loop time: 0.4

At the bottom of the dialog are 'OK' and 'Cancel' buttons.



6.7.5 Configuring the Transmitter and Pre-amplifier Pre-delays

The command `edscon` allows to set the pre-delays for the transmitters and preamplifiers. Currently, these are not used by the Fourier system. Instead, they are preset to reasonable values. These can be changed on the Fourier web page. This may be changed in future releases. For now, changing these values will not effect the hardware settings.



6.8 Setting the Magnet to the Right Field

The magnet is designed to be operated at 300.18 MHz with the current in the H0 Field Coil of the Shim System set to zero. At this frequency the Synthesizers are "cleanest" and least likely to create spurious problems.



When the lock is powered up from cold, the Field DAC is not automatically set.


In order to make sure the current in the Field Coil is set to zero

- Open the BSMS panel.
- Set the Field Value to 2048 and switch the lock to ON once.
- Switch it off again, while the magnet is bring charged to the right field and the cryo-shims are adjusted.

6.9 The Lock System

6.9.1 Working with the Lock System

The Fourier lock system is very similar to the BSMS lock. The BSMS Control Suite is fully functional and for normal operation, the differences are negligible.



As explained in the previous chapter, some of the parameters have a different value range. This applies mainly to the PID parameters. Increasing loop gain and loop time results in a more rigid coupling to the lock signal (fast regulation with a higher risk to see lock artifacts in the baseline next to intense resonance lines). Decreasing them results in less coupling (slow regulation, not suited to fast disturbances, but less baseline artifacts).

The *loop filter* is only applied on the lock data as they are displayed on the screen. The higher the value, the more data points are averaged to reduce noise.

6.9.2 Finding a Lock

If the lock is set up correctly, there is no need to search for the lock manually. The lock command will result in the lock system to read all parameters from the edlock file and to perform an automatic search and lock procedure.

During installation or in case the magnetic field changed by a considerable amount, it is necessary to find the lock signal for at least one solvent. This way, the field value defining the B0 shim current can be set correctly.

There are two possibilities to find the lock.

First approach: Sweeping the B0 field.





Figure 6.5: BSMS Panel

On the BSMS panel the sweep mode can be activated. The range (Ampl) and rate of the sweep can be adjusted. If there is no signal close to the current field value, the lock window only shows a noisy baseline. Use the **Field** value from the BSMS panel to search for the lock signal. Don't increase the value too fast because it is easy to miss the signal otherwise. The lock signal looks like this if the field value is correct: symmetric around the middle of the screen with both leading edges of the signal showing up. The latter must be achieved by adjusting the **Phase** of the lock.



Once, these parameters are set correctly, open the **edlock** window. Select **BSMS** in the menu and press Field Update. Leave the Phase set to -1, since this is stored by the lock control.

The screenshot shows the Edlock software window with the following table of solvent parameters:

Δ Solvent	Lock Power	Loop Gain	Loop Time	Loop Filter	Lock Phase
Acetone	-29	0.4	0.4	4	-1
C6D6	-20	0.4	0.4	4	-1
C6D6+Dioxane	-20	0.4	0.4	4	-1
CD2Cl2	-20	0.4	0.4	4	-1
CDCl3	-12	0.4	0.4	1000	-1
D2O	-20	0.4	0.4	4	-1
Dioxane	-20	0.4	0.4	4	-1
DMSO	-20	0.4	0.4	4	-1
EtOD	-20	0.4	0.4	4	-1
H2O+D2O	-12	0.4	0.4	1000	-1
MeOD	-12	0.4	0.4	4	-1
Tol	-28	0.4	0.4	1000	-1

At the bottom right of the window, there are two buttons: "Restore default" and "Close".

There is a problem with this approach which also applies to the BSMS. The lock sends a series of RF pulses with a rate of 1kHz to the 2H coil. This results in an excitation profile of a central peak and a number of additional peaks on a 1kHz grid. These slowly decay in amplitude as their distance to the central peak increases.

Setting the sweep to a large amplitude clearly shows the problem.



There is more than one *resonance* due to the excitation grid and it is not obvious which one is due to the central peak.

The rate of RF pulses on the BSMS is 6.66kHz. So, the effect is less obvious and the determination of the right sweep signal is much easier.

Second Approach: The "flock" AU Program

There is an alternate way of finding the correct lock field. Instead of using the lock instruction, it is possible to start an AU program called flock. flock asks for the lock solvent, creates a temporary 2H dataset and starts an acquisition to find the lock signal in the 2H spectrum. It compares the ppm shift of the 2H signal to the default value from the **edlock** table. The difference in ppm is then converted to a shift in the lock field. The new lock field is displayed in a pop-up window. Additionally, flock informs the operator how many solvent peaks have been found and if this number is in agreement with the number of resonance lines for the selected lock solvent. In case, the number of signals is different from the expected value, flock uses the peak with the highest intensity to calculate the lock shift. The new shift value is set in the Fourier and a normal lock procedure is started.






Figure 6.6: Result of Flock with one Solvent Peak



Result of flock if number of detected peaks is not compatible with number of expected solvent peaks. This can happen if a mixture of solvents is used or if the wrong solvent has been chosen. In any case, it leads to a field lock. The ppm scale may be shifted if the wrong solvent has been chosen, though.

Any time, the lock field has changed by a large amount, the following message is shown:



Pressing **OK** updates the BSMS field value in **edlock** table. A subsequent lock procedure should be able to find the lock signal in the usual way.

Fine-tuning the flock AU program.

In the AU program, one can find the line:

```
static const double scale = 19.79;
```

It defines the ratio between the ppm scale and the BSMS field values. In principle, any Fourier should be close to that value. It is recommended, though, to calibrate it for any instruments as the shim current boards are not completely identical to each other.


To determine this value, one has to acquire two proton spectra. For this the Lock should be OFF. After acquisition of the first spectrum, one changes the lock field by 1000 units (+/- doesn't matter). With the new field value another spectrum is acquired. Now, determine the ppm shift of the two spectra in respect to each other, e.g. in the dual display mode. This value should be close to 20ppm. Please correct the line of the AU program shown above to reflect the accurate value (take the absolute value, no negative sign).

Third approach: Acquisition of a proton spectrum


An easier and unambiguous approach is to acquire a proton signal using a large spectrum width SW. The easiest sample is *Doped Water* with just one large resonance line from H₂O. The solvent is D₂O. Therefore, the lock frequency in ppm on the 2H scale is very close to that of the proton spectrum on the 1H scale (= 4.7ppm). First try to lock on D₂O in order to transfer the lock parameters for D₂O to the lock system. The lock will fail because we are in a situation there the signal is outside of the lock window. Otherwise, we wouldn't need to search for the lock. With the Sweep set to off, change the field value such that the 1H resonance line is at 4.7ppm. This is also the right field value for the lock signal. Now, inspecting the lock window, it is possible to set the phase correctly. Save both values in the **edlock** window to finish the procedure.

The procedure is illustrated using Cyclosporine as a sample. There is always a resonance line from the protonated lock solvent because chemical exchange is easy between 1H and 2H. In cyclosporine, it is the largest signal on the left-handed side of the spectrum resulting from C6D₆H.

If the field value is off, it may look like this.




Here is a zoom to the region of interest. The resonance line is at 26.22 ppm.



From the edlock table, one can read that the C6D6 chemical shift for 2H is 7.16ppm. So, use the field value to shift the proton resonance line to that value on the 1H frequency scale. A rough value for the field shift is 20ppm/1000 field units. It depends a bit on the actual shim board but can be a good starting point for all Fourier systems.

There is no need to be exact at this point. A spectrum like this is close enough. Using the resulting field value, the Fourier must be able to lock the field in the automatic mode.



6.10 Converting the Parameter Set to the Actual BF1

Since the shimming process may have an effect on the B0 field, it may be necessary to redefine BF1 after Shimming and Lock have been optimized. The ideal BF1 is close to 300.18 when the Lock Field DAC is close to 2000 (± 500).

When you are satisfied that you have a good BF1 and the LOCK DAC is near its center, repeat **CF** and **expinstall** with the new BF1.

6.11 Tuning and Matching the Probe

6.11.1 General Information about the Fourier Probe

The Fourier system is delivered with a special probe. It is a modified DUL probe with 13C and 2H on the inner coil and 1H on the outer coil. The special feature is a coil design which is much less susceptible to the salt content of the sample. Tuning and matching are much less influenced by a change of the sample. Therefore, the probe is tuned and matched once and this calibration is not changed on a running system. The side effect is a lower S/N ratio by about 15% compared to a normal DUL probe.

There are no tuning and matching rods on the bottom of the probe. The tunable capacitors are accessible via small holes in the upper part of the probe body. Use the screwdriver coming with the probe to tune and match the probe before installing it into the shim system.

The wobb routine of the Avance system does not work with the Fourier. Instead, there is a special wobb program which must be started outside of TopSpin.

See also chapter [Additional Test Applications \[▶ 100\]](#) for details.

6.12 The Gradshim Command

In order to enable automatic shimming, the probe is equipped with a z gradient. Unfortunately, TopShim doesn't fit to the pulse programmer structure of the Fourier system. Therefore, Gradshim is used instead.

It should be noted that only z gradients can be used for shimming. The focus is not so much on the ease of installation but more on the automatic shimming when the spectrometer is running in automation.

Gradshim works best when the shim map was made with an already well shimmed magnet. You can also use a shim map from another system to get started, since the Fourier Shim Control and BOSS 1 Shim Systems are fairly consistent.

Also important for the success of gradshim are:

- The gradient echo is nicely centered in the spectral width - check by loading and running the experiment gradshim1d2h_f.






Figure 6.7: Gradient echo

- If this does not look like in the figure above, it could be, that the parameter sets have not been converted to the correct BF1.
- The correct 2H pulse length. This is normally around 50us at 40 W and can be checked by a zg experiment, where 2H is selected as the observe nucleus. If the lock is still running while you do this, there will be spikes in the 2H spectrum!

The command **gradshim** opens a new window.



If **gradshim** is done for the first time on a system, the profiles of the various shim gradients must be evaluated. (open Setup in the main window).



This has only to be done once. If the shim profiles are available, it is possible to *Start Gradient Shimming*.



IMPORTANT: To start gradshimming using the **shim** button of the flow bar and also from IconNMR, it is necessary to store the current configuration for automation. (Setup -> store for automation). The flow bar or ICON will call either a 1H or a 2H gradshim, depending on what was selected when the setup was stored for automation.

For further information on **gradshim** please check the respective manual.



Gradshim uses different parameters on a Fourier. Be sure to define the parameter sets gradshim1d1h_f or gradshim1d2h_f when creating a shim map.

6.13 Installation Acceptance



All the for the Fourier 300 necessary Hardware Tests (HWT) were carried out during console end testing and are not intended to be repeated at the customer site.

The following ATP tests are carried out on the FOURIER during customer acceptance

- Determination 90 degree 2H observe pulse
- Determination 90 degree 1H observe pulse
- Determination 90 degree 13C observe pulse
- Determination 90 degree 1H decoupling (CPD) pulse
- Determination 90 degree 1H high power decoupling pulse
- Determination 90 degree 13C low power decoupling pulse
- Line shape test for 1H with rotation
- Resolution test for 1H

- Sensitivity test for 1H
- Line shape Test for 13C with rotation
- Resolution test for 13C
- Sensitivity test for 13C (ASTM)
- Sensitivity test for 13C (EB)
- HSQC experiment with Z-gradient (console test)
- COSY experiment with Z-gradient (console test)
- DEPT-90 experiment
- DEPT-135 experiment

For the ATP Tests you will need the following samples

The three Test Samples which were shipped with the system.

- 3% CDCl₃ in Chloroform (line shape)
- 0.1% EB in CDCl₃ (1H SINO)
- ASTM (C13 SINO)

To shim the magnet and perform the ATPs test you will require to bring some additional samples:

- Urea in DMSO (Pulse Length adjustment)
- 100mg or 1g Cholac in CDCl₃ (DEPT)
- Cyclosporine in C₆D₆ (2D Examples)

There are only three hard specs for the Fourier required for a customer sign off:

Lineshape using 3% CDCl ₃ in Acetone	0.5/8/16 Hz (non-spinning)
1H Sensitivity (0.1% EB in CDCL ₃)	125:1 (after gradshim)
C13 Sensitivity (ASTM)	90:1 (after gradshim)

Table 6.1: Hard Specs for the Fourier

The remaining ATP Tests will implicitly demonstrate the system remaining functions:

- Water Suppression
- Sample Spin & Lift
- Auto-Locking on various solvents
- Gradient Functionality (2Ds & Gradshim)
- Gradshim reproducibility
- Lock and Temperature Stability (COSY & HSCQ)
- * Baseline & Spectral Purity (DEPT)

* Certain artifacts are normal and have been considered acceptable for the intended use of the Fourier:

- A moderate amount of T1 noise in 2D experiments
- Some "digital noise" on strong signals that require RG < 10. Use zg30 instead zg of if that happens, in order to be able to increase RG
- Some vibration noise close to the feet of very narrow lines (e.g. lineshape samples). This can be improved with vibration damping if required.

- * Some very small spurious spikes in the first few scans of C13 spectra. These should average out after NS > 64
- *Asymmetric smileys* in some 1D spectra with high intensity. This can also be seen as reduction in the level of noise about 200 Hz from the right edge of some 1D and 2D spectra. If this part contains spectral information, simply increase the spectral width.

* If you obtain reproducible "spurious spikes" in 1H spectra it may be due to an unfavorable BF1 frequency. The Fourier DDS are "clean" for at least 200ppm round 300.18 MHz

6.14 Main Differences between Fourier & Avance Systems

- The Lock works at 1 kHz instead of 6 kHz. That means that sidebands are closer together. Also the lock parameters have different units to BSMS Values.

P:	0.1 to 0.8 depending on solvent. Higher values mean more proportional gain.
I:	0.1 to 0.8 depending on solvent. Higher values mean more integral gain.
Power:	-12 to -28 dBm depending on solvent. Note: The values shown in the BSMS Panel are attenuations from 0dBm.
Filter:	4 for strong solvents, 1000 for weak solvents. The numbers are not the Bandwidth in Hz but the number of averages per second. 1000 is equivalent to 1 Hz, 4 is equivalent to 250 Hz.
Field:	0 to 4095. 2048 means no current flowing in H0 Coil. Ideal value is between 1500 and 2500, to allow for magnet drift. The BSMS field has the opposite polarity, i.e. increasing the value of the field DAC by 1000 steps will cause a shift left (i.e. increase the frequency) by 20 ppm.

Table 6.2: Typical Values for the Fourier

- Shim Values go from -32768 to 32768. If converting from a AVANCE shim file, you will need to halve the values to get a good starting point.
- No PID Self Tune for Lock, VTU and Spin Parameters. Since the Fourier is only sold with one probe, the defaults if operated at 4 bar are usually good enough for all experiments.
- Different Pulse Programmer Language
- Digitizer Type must be set to "DRMS"
- Aq.mod must be set to qsim
- Fn_Mod must be set to TPPI or QF
- BC-Mod must always be set to OFF
- GPNAM must always be RECT.1

7 Configuration

Each of the three Ethernet Devices has an ARM CPU with a web server for service, BIS setup and diagnostic functions.

In order to prevent accidental damage, the web server access is password protected:

User: **root**

Password: **BRUKER (in capitals !!)**

You can connect to the web server using any browser, but Firefox is recommended in case you want to save the password.

7.1 Network Settings

The arm CPU receives its IP address by different methods, selected by jumpers on the dip switch shown below. The default for use with TopSpin is DHCP, i.e. (X, ON, OFF, OFF (X = don't care)

If for any reason you cannot connect to a device, it may be necessary to set the jumpers to a default address (e.g. X,OFF,OFF,ON for 192.168.1.1)



Both the lock board and the main board are based on the minispec architecture, so if they are set to X,OFF,OFF,OFF will both default to 192.168.1.200. The VTUs default is 192.168.1.109.

It is also possible to assign a fixed user specific IP address.

Main -> support -> network configuration

Parameter	Software	EEPROM	Modify EEPROM
IP address:	192.168.1.1	149.236.42.97	<input type="text" value="149.236.42.97"/>
Subnet Mask:	255.255.255.0	255.255.255.0	<input type="text" value="255.255.255.0"/>
Standard Gateway:	0.0.0.0	0.0.0.0	<input type="text" value="0.0.0.0"/>
MAC Address:	0000AD08FF12	0000AD119C12	---
SW2 jumper settings:	OFF, OFF, ON	ON, ON, ON	---

Figure 7.1: Fixed User Specific IP Address


This address will be used if the jumpers are set to X,ON, ON, ON.

The jumper settings are described on each web server on

Main -> help-> power up jumper

Fourier

Power up jumper



2 / 6	3 / 7	4 / 8	DHCP	BOOTP	MAC address	IP address
						no hardware found
						192.168.1.100
						Hall
						192.168.1.104
						Signal Channel
						192.168.1.106
						Bridge
						192.168.1.107
						ENDOR
						192.168.1.108
OFF	OFF	OFF	no	no	defined by serial number	VTU
						192.168.1.109
						Teslameter
						192.168.1.110
						Signal Channel 2
						192.168.1.111
						Gradient Unit BMC20
						192.168.1.114
						SPU
						192.168.1.117
						Cryo Power Supply
						192.168.1.150
						MiniSpec
						192.168.1.200
OFF	OFF	ON	no	no	0x0000AD0BFF12	192.168.1.1
ON	ON	OFF	yes	yes	defined by serial number	from DHCP/BOOTP server
OFF	ON	OFF	no	yes	defined by serial number	BOOTP server
ON	OFF	OFF	yes	no	defined by serial number	DHCP server
ON	ON	ON	no	no	defined by serial number	user defined

Figure 7.2: Power Up Jumper

7.1.1 Resetting the CPU

Any change in network configuration, firmware of FPA content are only valid after rebooting the CPU. There are three ways to reset the CPU:

By

- Powering up
- Pressing the **reset** button on the CPU
- Pressing the **Hardware reset** button in the web server on

Main -> support -> network configuration


This button also appears after a successful firmware update.

7.2 Firmware Update

The firmware of all the intelligent devices can be upgraded during run-time via the web-browser. The mechanism for the Ethernet Devices differs to that of the CAN-Bus Devices.

7.2.1 Ethernet Devices (Main Board, Lock Board, VTU)

Before changing Firmware the current firmware version must be identified. In Ethernet Devices this is always shown at the bottom of each web server page.



To download a new firmware select


Main -> support -> CPU software update

1. Select the firmware to be uploaded using the **browse** button (in some versions called "Durchsuchen").
2. Press **Upload** to upload firmware to CPU.



This takes a few seconds and is finished when the message appears: **File ... uploaded successfully!**

3. Press **Flash** to write the uploaded file into the Flash memory of the CPU.




This takes about a minute and is finished when the message appears: **Successfully Programmed Flash.**

4. Reset the CPU by pressing the **Hardware Reset** button.



5. Confirm the new version after the CPU has been reset.



7.2.2 CAN-Bus Devices



The procedure is the same for all devices.


The procedure is the same for all devices, however

- since the **shim control board** and the **GAB/2 Board** are slaves of the fourier main board, the firmware update of these devices is initiated form the **main board web server**.
- since the **peripheral control board** is a slave of the VTU, the firmware update of this devices is initiated form the **VTU web server**.

To start the procedure select

Main -> support -> AVR Controller software update

1. **Select Device** and Board Number and press **SET**.



GAB/2 Backplane and Peripheral Control the Board Number is always “Board 1” Since the shim control comprises two boards, this procedure has to be carried out twice. Once for Board1 and once for Board2.

2. Select the EEPROM file (file ending *.eep) to be uploaded by clicking on **Upload EEPROM Data**.

VTU
Upload EEPROM Data

AVR Controller:

AVR AT90can128	
EEPROM Memory Size	4 kBytes

Write generic software:

*.eep File :

This procedure takes less than a second and is finished when the message appears: **Flash EEP File**.

3. Press **Flash Eep-File** to write the file to the EEPROM.

VTU
AVR Flash EEPROM Data

File: **peripheral.eep** of size: **103** has been uploaded successfully !

This procedure takes about 30 seconds and is finished when the programming menu reappears.

VTU
AVR Software Update

Application:

Application Code	PERIPHERAL	Board Number	Board 1
------------------	------------	--------------	---------

Software Update:

- [Upload EEPROM Data](#)
- [Upload Flash HEX - File](#)

Hex File

4. Repeat this procedure (Upload & then flash) for the Hex file.

VTU
AVR Flash Hex-File

File: **peripheral.hex** of size: **69236** has been uploaded successfully !

This procedure takes a few minutes and is finished when the following message appears:

VTU
AVR Software configuration

Download code is 0: Successfully programmed flash.
You can reset the hardware or flash additional software now.

5. **Reset** the Master CPU (in this example the VTU) by pressing the **Hardware Reset** Button.

Shim Control: Main Board: main -> shim -> base function

6. **Check** that the correct version has been successfully updated.

Version Board	2.9	Read
Version Board	2.9	Read

GAB/ 2: Main Board: main -> fourier --> GAB/2 System -> GAB/2 Menu -> GAB/2 Control

GAB/2 Status		
Relay Mode		on
Analog Path		on
GCON Parity Check		inactive
GAB/2 Backplane Firmware Version		V00.04
GAB/2 Backplane FPGA Version		V0.01

[Main](#) | [MiniSpec](#) | [GAB/2 System](#) | [GAB/2 Menu](#) | [GAB/2 Control](#)
[Main](#) | [Support](#) | [MiniSpec](#) | [Shim](#) | [Help](#)

Peripheral Control: VTU: main -> VTU -> peripheral control -> base functions

Messages:

Function	Parameter	Method
BeAlive	Be alive	Get
Version	0.016	Get
LED	on/off	Set

7.3 FPGA Update

The FPGA Content of the Ethernet devices can be upgraded during run-time via the web-browser. The main board requires an external tool to stream the FPGA content, since the memory of the embedded web server CPU is not large enough to hold an entire FPGA image.

7.3.1 Lock Control and VTU

With these two devices, the FPGA is updated using the web server.



Before a start the current FPGA version must be correctly identified.

1. Open **main->device->base functions** and read out the address of register 0x04.

Lock Board
MiniSpec Base Functions

Write Command: Read Command:

Register	Value
0x00	0x00

Register	Value
0x04	88 = 0x00000058

2. Begin the update procedure using **main -> support -> XILINX CPLD/ FPGA Update**.
3. Select the corresponding FPGA image (file ending *.xsvf) and press **Upload Xilinx XSVF-File**.

Lock Board
Upload Xilinx XSVF-File

File: browse

[Main](#) | [Support](#) | [Upload Xilinx XSVF-File](#)
[Main](#) | [Support](#) | [MiniSpec](#) | [Help](#)

Version: M4.69 - Build: 1490
 CPU (V2) SN: 1947 - HW SN: ---
 ©BRUKER BIOSPIN GMBH

This takes a few seconds and is finished when the **Start** window appears. 4. Press **Start**. This procedure takes about **25 minutes !!!**

Xilinx JTAG programming

File: **minispec_v84.xsvf** of Size: **2719250** has been uploaded successfully !

Disconnect the JTAG download cable !

It is finished when this message appears.

5. **Reset** the CPU by pressing the **Hardware Reset** button in the menu

Xilinx programming result

Download code is 0: Successfully programmed device.

6. **Verify** that the new version has been correctly loaded by reading register 0x04 (see above).


7.3.2 Main Board

The main board requires an external tool to stream the FPGA content, since the memory of the embedded web server CPU is not large enough to hold an entire FPGA image. The procedure is therefore different (and simpler) but does need a special software tool called **Allegro Flash Tool (Flash.exe)**.



Before a start the current FPGA version must be correctly identified. In the case of the main board this is seen at the bottom of the main web server page.

1. Identify the current FPGA version.




2. Start Allegro **Flash Tool** (Flash.exe), enter the IP address of the main board, and select the FPGA image (file ending *.lz)

3. Press **Flash** to start the procedure.

This procedure takes about **35 minutes !!!**



It is finished when the following message appears:



4. **Reset** the CPU by pressing the **Hardware Reset** Button in the menu **Main -> support -> network configuration**.

5. **Verify** that the new version has been correctly (bottom of main web page).

7.4 Device Specific Settings & Tools

Apart for the generic functions for network configuration, Firmware & FPGA Update each web server has some individual device specific functions such as:

- BIS EEPROM
- Device Identification
- Device Setup Parameters
- Diagnostic Functions

7.4.1 Main Board

7.4.1.1 BIS EEPROM Settings

To recognize the main board as Fourier component, the following entries need to be made in the field **NAME group** part of the BIS EEPROM. This is accessed via:

Main -> minispec -> BIS information

PRODUCTION group:

Product ID:	E1451110
Serial number:	NE9999
ECL:	0.0
Variant:	
Location:	
Date:	14/04/2010

Apply

NAME group:
Enter device identifier ("MQD"; DreamSpec).

Unk\IMQD

Apply

INITIALIZATION group:

Mode (MINISPEC or IMAGING):	MINISPEC
Internal use 1:	
Internal use 2:	
Internal use 3:	
Internal use 4:	

Apply



The names *minispec* and *dreamspec* are historical and are used internally only.

Next, further down on the same page, the **frequency range** has to be defined. For 300 MHz this is:

Apply

DreamSpec settings:

Operating Frequency Min:	295
Operating Frequency Max:	305
NMR Separation Frequency:	35
DDS Sync Delay Mod1:	0x28000050
DDS Sync Delay Mod2:	0x28000048
DDS Sync Delay LO:	0x28000050

Apply



The three **DDS Sync Delay** parameters are board specific delay line values that are determined to make sure that all three DDSs are synchronized and should not be changed!

Next the dependencies and Polarities of the **service pulses** have to be defined. There are four programmable service pulses, and they are used to gate transmitters and preamplifiers. For each pulse the HF pulse that it should be connected to, the polarity and a pre- and/or post-delay (in multiples of 40 ns) can be defined. The default settings are:

DreamSpec transmitter pulse (service pulse 0):

Transmitter Master (MOD1 or MOD2):	MOD2
Transmitter Blank Polarity, ST2 (0 = active high, 1 = active low):	0
Transmitter Blank Predelay / 40ns:	25
Transmitter Blank Postdelay / 40ns:	0

1H Transmitter gate
(with 1 us predelay)

Apply

DreamSpec preamp pulse (service pulse 1):

Preamp Master (MOD1 or MOD2):	MOD2
Preamp Protection Polarity, ST8 (0 = active high, 1 = active low):	1
Preamp Protection Predelay / 40ns:	0
Preamp Protection Postdelay / 40ns:	0

1H Preamp Gate
(inverted)

Apply

DreamSpec Mod2 pulse (service pulse 2):

Mod2 Master (MOD1 or MOD2):	MOD1
Mod2 Gate Polarity, ST4 (0 = active high, 1 = active low):	0
Mod2 Gate Predelay / 40ns:	25
Mod2 Gate Postdelay / 40ns:	0

C13 Transmitter gate
(with 1 us predelay)

Apply

DreamSpec SP3 pulse (service pulse 3):

SP3 Master (MOD1 or MOD2):	MOD1
SP3 Gate Polarity, ST9 (0 = active high, 1 = active low):	0
SP3 Gate Predelay / 40ns:	0
SP3 Gate Postdelay / 40ns:	0

C13 Preamp Gate

Apply



The preamp gates are not protection pulses. The preamps are self protecting using diodes, and the gate is required in order to ensure that low power pulses can pass the protection diodes.

7.4.1.2 Board Temperature Control

The temperature sensitive parts (Master Oscillator, Receiver Gain Stages, Filters, etc.) are all within a temperature controlled zone on the main board, which is thermally isolated with a blue foam cover.

The temperature is controlled by a digital PID algorithm which runs in the background on the embedded CPU. The PID settings can be adjusted using:

Main-> minispec-> electronics temperature control

Fourier

Electronics temperature control

AD controller 2 (NTC):

Chn 1	Chn 2	Chn 1 / Chn 2
6881640	6879890	1.0003

actual temperature

Read

Heater parameter:

Heater PID	Set value	Parameter
Proportional band	3	10000
Integral time	1483	20
Derivative time	0	0

Apply

temperature setpoint heater power

Temperature setpoint	Heater register	Heater power [%]	LED (counter)
1.0000	1485	36.3	TBD

Apply

A temperature setpoint of 1.0 is equivalent to about 40°C. It should be adjusted so that the heater power is between 30 and 50 % when stable conditions are reached. Adjust the setpoint if necessary (lower value -> higher temperature. The setpoint is raised by 1°C when the setpoint value is lowered by 0.05, and vice versa).

PID Values should not need to be changed. The values are stored in the EEPROM and are valid after power up.

7.4.1.3 On-Board Voltages

The main board receives only 28Vdc from the central 28V PSU and generates all the other voltages it requires on board. You can test if all the voltages are correct under:

Main-> minispec-> AD Controller

20	1BD	37.255	[1BD]	[1BD]
21	28V_PA2	26.754	28	V
22	28V_PA1	22.751	28	V
23	4V	4.017	4	V
24	2V5_FPGA	2.417	2.5	V
25	3V3_FPGA	3.413	3.3	V
26	3V3_DDS	3.383	3.3	V
27	3V3_ADC	3.294	3.3	V
28	5V_DIGITAL	5.177	5.2	V
29	5V_ANALOG	5.141	5.2	V
30	28V_EXTERN	27.610	28	V
...	---	-----	-----	-----

7.4.1.4 Pulse Programmer Memory

Unlike the IPSO, the Fourier keeps its entire Pulse Program (incl. decoupling sequences, loops, 2D evolutions) in memory at once. Since the memory is limited, some pulse programs may exceed the memory capacity. In order to check if this is the case, you can check the memory usage using:

Main-> minispec-> memory

Pulse Programmer	Command	Table
Modulator 1	10.64 %	3.22 %
Modulator 2	1.56 %	3.22 %
Receiver	1.56 %	3.17 %
Gradients	0.00 %	0.00 %
Service channel 1	1.56 %	3.13 %
Service channel 2	1.56 %	3.13 %
Service channel 3	10.64 %	3.13 %
Service channel 4	10.64 %	3.13 %

Read

7.4.2 Shim Control

Since shim control is a slave of the main board, the shim specific web server functions are found on the main board web server under:


Main-> shim



The Fourier has 20 coils for shimming. One Shim Control Board can handle 12 coils. Therefore we need two Shim Control Boards. The first Board is represented by the Menu item “Coil 1”, the second Shim Control Board is represented by the Menu item “Coil 2”.

Base Functions

Here you can see the Firmware version of the two boards. If you want to test the Can-Bus communication you can toggle the 4 LED on board. If you mark the checkbox you can activate the LED. Below the item *Temperature* you can see the temperature of the two boards.



Coil 1

On this page you can see the 12 coils which are controlled by the first Shim Control Board. The resolution of the DAC which set the current to the coils is 16 bits. 0x000 is the maximal negative current and 0xFFFF is the maximal positive current. This means the half 0x8000 is zero current. In the column “current” you can read the current through the coil. The next column ratio shows you the ratio of the current which flows the coils and the current which should flow. This ratio should be near 1.



If the current is close to zero, sometimes an error *insufficient current* is shown. If you are not sure if a shim is responding, set a larger current (> 100mA) temporarily and check that the error is fixed.

Coil:	name	DAC		Current [mA]	Ratio
X		0x8000	Set	-0.39	1.00
Y		0x8000	Set	-0.42	1.00
Z		0x8000	Set	-0.16	1.00
XY		0x8000	Set	-0.52	1.00
XZ		0x8000	Set	-0.42	1.00
YZ		0x8000	Set	0.00	1.00
R ²		0x8000	Set	-0.15	1.00
x ² - y ²		0x8000	Set	-0.40	1.00
x ²		0x8000	Set	-0.23	1.00
y ²		0x8000	Set	0.17	1.00
z ²		0x8000	Set	-0.16	1.00
R ⁴		0x8000	Set	-0.13	1.00
H ⁰				0.00	

Write DAC-Values to FRAM
 Read DAC-Values from FRAM
 Get Dac values
 Set DAC-Values to coils

Apply

Write DAC-Values to FRAM:

The DAC-Values will be saved to the onboard storage (non-volatile RAM)

Read DAC-Values from FRAM:

If you want to get the DAC-Values you have saved in the FRAM, you can choose this option. Please remember, the values will only catch from the FRAM to the controller on the Shim Control Board. After this operation you have to do the next operation **Get Dac-Values**

Get Dac values:

This function fetches the active DAC-Values from the Shim Control Board controller to the Web page.

Set DAC-Values to coils:

This function sets the active DAC-Values in the Shim Control Board controller to all coils.

Apply:

This button will apply the function you have chosen with the checkboxes.

Coil 2

On this page you can see the 8 coils which are controlled by the second Shim Control Board. The resolution of the DAC which set the current to the coils is 16 bits.


Coil:	name	DAC		Current [mA]	Ratio
	Z5	<input type="text" value="0x8000"/>	<input type="button" value="Set"/>	-0.12	1.00
	YZ2	<input type="text" value="0x8000"/>	<input type="button" value="Set"/>	-0.32	1.00
	YZ3	<input type="text" value="0x8000"/>	<input type="button" value="Set"/>	-0.21	1.00
	Z6	<input type="text" value="0x8000"/>	<input type="button" value="Set"/>	-0.13	1.00
	(X2-Y2)Z	<input type="text" value="0x8000"/>	<input type="button" value="Set"/>	-0.23	1.00
	XZ3	<input type="text" value="0x8000"/>	<input type="button" value="Set"/>	-0.37	1.00
	XYZ	<input type="text" value="0x8000"/>	<input type="button" value="Set"/>	0.04	1.00
	XZ2	<input type="text" value="0x8000"/>	<input type="button" value="Set"/>	-0.04	1.00

Write DAC-Values to FRAM
 Read DAC-Values from FRAM
 Get Dac values
 Set DAC-Values to coils

7.4.3 GAB/2 Interface


Since the GAB/ Backplane is a slave of the main board, the GAB specific web server functions are found on the main board web server by selecting:

Main-> minispec -> GAB/2 System



In order for the GAB/2 Backplane main board to be recognized as a Fourier component, the following entries need to be made in NAME group part of the BIS EEPROM. This is accessed via:

Main-> minispec -> GAB/2 System -> BIS Information



You can check the presence of the GAB/2 Backplane and the GAB/2 Gradient amplifier using:

Main-> minispec -> GAB/2 System -> Base functions

GAB System Configuration:

GAB System	Device Status	Connection Status	Operation
GAB System Backplane	present	connected	<input type="button" value="Check"/>
GAB/2 Board	present	connected	<input type="button" value="Check"/>

You can check if a LVDS link is present using:

Main-> minispec -> GAB/2 System -> GAB/2 Menu -> GAB/2 Control

GAB/2 Command	
GAB Mode	<input type="button" value="ON"/> <input type="button" value="OFF"/>
pending (waiting for LVDS link)	
GAB/2 Status	
Relay Mode	off
Analog Path	off
GCON Parity Check	inactive
GAB/2 Backplane Firmware Version	V00.04
GAB/2 Backplane FPGA Version	V0.01

7.4.4 Lock Board

In order for the lock main board to be recognized as a Fourier component, the following entries need to be made in NAME group part of the BIS EEPROM. Select

Main -> minispec -> BIS Settings

NAME group:

Enter device identifier ("MQ1": mqOne).

Unit:

Select the main information and setting page by

Main -> minispec -> Lock

This page is divided in 4 parts:

- properties for the pulse programmer
- the current status of the lock signal
- Fine DAC search function
- Coarse DAC sweep function

The dynamic parameters (lock power, frequency, phase, loop gain etc.) of the lock panel are set by the **edlock** table. The values circled in red are base settings and need to be set correctly for proper lock operation.

An explanation of the parameters and functions follows:

Configuration

Pulse properties:

pulse length (ms)	dead time (ms)	number of points	recycle delay (ms)
0.100000	0.020000	4000	1.000000
NMR frequency (MHz)	Gain (dB)	Attenuator (dB)	Magnitude Lock Lost
46.079787	80	20	100000
P-gain	I-gain	Integral Setpoint	Freq. Calibration (Hz)
0.100000	0.100000	0	125
Display Scale Factor	Field DAC Small Field DAC Wide	moving average	phase angle(deg)
2.17	647 1508	4	26.0
<input type="checkbox"/> intern = unchecked extern = checked <input type="checkbox"/> Self Test <input type="button" value="Start/Update"/> <input type="button" value="Stop"/>	<input checked="" type="checkbox"/> auto range <input type="button" value="P-Control ON"/> <input type="button" value="I-Control ON"/> <input type="button" value="Control Off"/>	<input type="button" value="operate +P +I"/>	<input type="button" value="Phase Correction"/>

DAC Value	Magnitude	real	imaginary
602	23319290	4829	7
<input type="button" value="Show"/>			

Field DAC center	Sweep width[steps]	Step size	Dwell time[ms]
2047	4000	40	500
Magnitude Max	DAC Val	Minimum Magnitude	
0	2047	100000	
<input type="button" value="Search"/>			

Pulse Length:

This value describes the time for pulses. **Default: 100 us**

Dead Time:

The dead time is the waiting time after the pulse before sampling begins. **Default: 20 us**

Recycle Delay:

This value represents the repetition rate for the pulse programmer. **Default: 1ms** (i.e. 1 kHz repetition)

Number of Points:

The number of data points averaged by the FPGA within the recycle delay. The maximum value is $(\text{recycle delay} - \text{pulse length} - \text{dead time}) / 12.5 \text{ MHz}$. **Default: 4000 points.**

NMR Frequency:

This is the 2H NMR frequency including chemical shift offset. There is no default as this value is set by **edlock**.

Gain:

This is the gain for the receiver path. 40 dB means no gain in relation to the pre-amplifier. The range is from 40dB to 106 dB. **Note** that this is NOT the gain set by the lock display panel which is a pure digital display factor. The hardware gain remains constant and should be set such that weak solvents have sufficient signal to reach the minimum lock and search levels, and low enough not to clip or change phase with strong solvents. **Default: 80 dB.**

Attenuator:

This adjusts the lock power level. At full power (0dB attenuation) the lock power is approximately 1mW (0dBm). The attenuator has a range 30dB in 1db steps. There is no default as this value is set by **edlock**.

Magnitude Lock Lost:

(Default 100000)

If the Magnitude of the lock signal (as seen in "show") falls below this limit, the **LOCK LOST** button in the TOPSPIN BSMS panel will go red show a weak signal or indicate that the lock has been lost.

P-Gain:

This is the proportional gain of the digital PI control. Higher Values mean tighter & faster control but also more lock noise and chance of oscillations. The lock will be stable up to values of around 0.8, but is normally set to around 0.1 by the **edlock** table

I-Gain:

This is the integral gain of the digital PI control. Higher Values mean tighter & faster control but also more lock noise and chance of oscillations. The lock will be stable up to values of around 0.8, but is normally set to around 0.1 by the **edlock** table

Integral Setpoint:

This is setpoint for the imaginary term during lock and should always be set to 0.

Frequency Calibration:

This value can be used to calibrate the Lock. This is done using a known sample. Increasing this value by 1 step will shift the 1H spectrum by approx. 6Hz to the left. **Default Value: 125**

Display Scale Factor:

This is a floating point number that scales the display gain. A value of 1.0 is equivalent to 60 dB in the BSMS panel.

Field DAC Small:

The Lock is controlled by two DACs. A 12-bit coarse DAC is used for range selection and a 12-bit fine DAC which does the actual field regulation. The resolution of the fine DAC is about 0.025 Hz per DAC Step.

A DAC value 0 means maximal negative current and 4095 means maximal positive current. 2047 means zero current. When the control loop is closed, this value is under control of the lock software.

Field DAC Wide:

The resolution of the wide/coarse DAC is about 50 Hz per DAC Step. When auto ranging is ON (default), this value is under control of the lock software. In case of magnet drift, when the fine DAC gets to below 500 or above 3500, the coarse DAC is incremented or decremented by one step, to bring the fine DAC back nearer the center of its operating region

Moving Average:

This is a moving average of lock monitor signal. The range is from 1 to 1000. There is no default as this value is set by **edlock**. A value of 1 represents a bandwidth of 1000 Hz, and a value of 1000 represents a bandwidth of 1 Hz.



Note that this is used only for the display of the lock monitor and not for the regulation. The regulation bandwidth is only affected by the P & I gains.

Phase Angle:

This is the Phase correction for the real and imaginary part of the lock signal. The Lock Level display in the front end TopSpin displays only the real part of the signal. The imaginary part is used to regulate the field. There is no default as this value is set by **edlock**.

Phase Correction:

In order to be able to lock, the real part needs to be maximum and positive and the imaginary needs to be close to zero. Pressing this button automatically calculates the phase required for this condition.

Start / Update:

Starts the pulse programmers and sets all parameters to the hardware. This is the OPERATE mode, which generates lock signal, but not yet close the control loop.

Stop:

This stops the pulse programmer and data collection. The lock is now in STANDBY mode and no RF Pulses are generated.

Auto Range:

In case of magnet drift, when the fine DAC gets to below 500 or above 3500, the coarse DAC is incremented or decremented by one step, to bring the fine DAC back nearer the center of its operating region. **Default: ON**

P-Control ON:

Enables the proportional part of the control loop (-> Operate +P). The fine DAC value will be the proportional gain times the difference between the imaginary part of the signal and zero.

I-Control ON:

Enables the Integral part of the control loop (-> Operate +I). The fine DAC value will be the integral gain times the accumulated difference between the imaginary part of the signal and zero.

When the Lock is enabled from within TopSpin it always enables BOTH the proportional and integral gain (-> Operate +P +I)

Control Off:

Switches off the control loop, but continues pulsing and data acquisition (-> Operate).

Standby/Operate:

This part of the Web page is a status field.

Status	Description
Standby	The Lock board is in Standby mode. No pulses are being generated, no data acquisition.
Operate	The pulse programmer is pulsing and data acquisition is active working. The PI control is OFF
+P	The PI control works with the P-Gain
+I	The PI control works with the I-Gain
Sweep on	The device is in the Sweep mode
Self test intern	The device is in the internal Self Test mode
Self test extern	The device is in the external Self Test mode

Table 7.1: Status and Description

Self Test:

The can self test the lock board by rerouting its output pulse to its receiver input. This creates a lock monitor signal of constant amplitude and almost no noise, which responds to phase and receiver gain like a real NMR signal. You may need to reduce the receiver gain, if the signal is too large.

Internal / External:

There are two modes for the self test. In the internal mode, the modulator output is routed to the receiver input internally on-board. This can be done purely by software. In the external mode you have to physically connect the output to the input, but this has the advantage that the pulse attenuator can also be tested.

Show

Pressing show refreshes the current values of fine DAC, magnitude, real and imaginary values of the Lock signal.

DAC-Value:

Shows the current value of the Fine DAC. If the PI control is active, this is an indication of what the lock is actually doing. When the lock is stable, this value should be close to the center of its operating region (i.e. around 2000) and can vary by up to ± 10 counts (each count represents 0.025 Hz). When Auto-ranging is enabled, when the fine DAC gets to below 500 or above 3500, the coarse DAC is incremented or decremented by one step, to bring the fine DAC back nearer the center of its operating region

Magnitude

This is a "magnitude" value of the lock signal and corresponds to the sum of the squares of the real and the imaginary part. Since there is no square root, taken at this point, it is strictly speaking a power value.

Real

This is the real part of the lock signal, which is used for the lock monitor

Imaginary

This is the imaginary real part of the lock signal, which is used for the PI control

Search

This performs a search over the entire Fine DAC range looking for the maximum signal in the 100 Hz or so (1H) covered by this range.

The Wide DAC will stay constant during this search function. It searches the range given by **sweep width (default 4000)** around the **Field DAC center (default 2047)** value in **steps (default 40)** number of steps, waiting dwell time (**default 500ms**) at each step.

At the end of the search (which with the above defaults takes about 20 seconds), the maximum signal found is displayed in **Magnitude Max** and the DAC value shown in **DAC val**. A search is considered successful, if the maximum value is above the value set in **Minimum Amplitude (Default 100000)**, and if the corresponding DAC value was between 500 and 3500 (i.e. not too close to its limits). If the search was unsuccessful, the **AUTO LOCK** procedure of TopSpin varies the coarse DAC in steps of ± 10 , ± 20 , ± 30 , ± 40 and ± 50 in order to find a good lock. If there still is no good lock result, the automatic routine flails, and the lock has to be found manually using the sweep mode.

Sweep

This function sweeps the wide DAC with a range of **sweep width (default 1000, i.e. about 5 kHz range)** around the currently set **Field DAC Wide** value, in steps of **step size (default 1)**.




Figure 7.3: Wide Coarse too high




Figure 7.4: Wide DAC too low




Figure 7.5: Wide DAC OK

The sweep mode can also be used to check and adjust Lock Phase.




Figure 7.6: Phase wrong




Figure 7.7: Phase Correct

7.4.5 VTU

In order for VTU the lock main board to be recognized as a Fourier component, the following entries need to be made in NAME group part of the BIS EEPROM. This is accessed via:

Main -> VTU -> BIS Settings

NAME group:
Enter device identifier ("vtu" for VTU).

Unit:

Select the main web page in **Main -> VTU -> Controller Parameters**.

VTU

VTU controller parameters

PID Controller:

Parameter description	Set value	Parameter
Proportional band	0.000000	0.500000
Integral time	0.000000	0.100000
Derivative time	0.000000	0.000000

Apply Changes

Temperature:

Setpoint [K]	Actual Temp. [K]
300.0	299.49

Apply Changes

Heater:

Heater [%]	Max [%]
0.000	100.0

Apply Changes

Temperature control:

Start Stop

This can be used to control the VTU independent from TopSpin, or to monitor its function. PID Parameters cannot be set here, this has to be done on a separate page (see below)

7.4.5.1 Changing PID Parameters

Select Main -> VTU -> PID Parameter

VTU

Upload PID parameters

Current parameter table

Temperature From [K]	Temperature To [K]	Gas Flow [%]	P term	I term	D term
250.00	350.00	0.00	0.50	0.10	0.00
0.00	0.00	0.00	0.00	0.00	0.00
0.00	0.00	0.00	0.00	0.00	0.00
0.00	0.00	0.00	0.00	0.00	0.00

The VTU can store different PID setting for different temperature ranges. Since the Fourier Probe VLC only has a narrow temperature range usually one range is sufficient.



The value for **gas flow** in this table must always be **0%** since this is only used in connection with an N2 Evaporator system.

Default Values are P=1.0, I=0.07 for Gas flows of 270 to 400 l/h.


The values can either be changed via the self tune function in TopSpin (not yet implemented in V 1.0) or uploaded manually. Since the firmware always expects a whole value table, you have so upload a file, even if you only want to change one value.

1. Create a text file with the following format:


250,350,0.0,1.0,0.07,0.0

(in this example the temp range in °C and 1.0 and 0.07 are P and I respectively)

2. Select this file with the browser and press upload.



3. Press start to save the values to the EEPROM of the VTU.



This only takes a second and is completed when the following message appears:




Successfully updated PID parameter table.

7.4.5.2 Setting the Heater Protection Threshold

In case the gas supply to the probe is missing, the heater switches off when the overheat sensor reaches the threshold safety limit. This can be set under:

Main -> VTU -> Thermocouple



7.4.5.3 Calibrating the Thermocouple

A 2 point correction (e.g. Ice water & Liquid N2) of the thermocouple reading can be done under:

Main -> VTU -> Thermocouple

Calibration:

Number	Actual Temperature [K]	Measured Temperature (WinAcquisit) [K]
1	<input type="text" value="77.00"/>	<input type="text" value="77.00"/>
2	<input type="text" value="273.30"/>	<input type="text" value="273.30"/>

7.4.6 Peripheral Control

Since peripheral control is a slave of the VTU, the specific web server functions are found on the VTU main board web server under:

Main-> VTU -> Peripheral Control




HARWARE CONTROL:

- [Base function](#)
- [BIS information](#)
- [Lift Control](#)
- [Spin](#)
- [Gas Flow](#)
- [Helium Level Sensor](#)
- [Board and Shim temperature](#)

7.4.6.1 Setting the Lift Parameters

When the sample lift is activated, the lift air valve is slowly opened from **start%** to **stop%** in steps of **step**, pausing **time** x 32ms at each step. Typical default values for an input pressure of 4 bar are shown in the example below. You can slow down the lift by increasing the time per step or reducing the step size. In order to save these values to the non volatile memory of the board, press **save to NVRAM**.

The lift can be tested from the web pageside using the up/down toggle function



7.4.6.2 Setting the Spin Parameters

When the sample spin is activated, the spin air valve is slowly opened with **offset%** and holds this value for a time given by the number of ticks (in the example below $3 \times 30 \times 32 \text{ ms}$ = approx. 3 seconds) before beginning the PI control. Typical default values for an input pressure of 4 bar are shown in the example below.

The starting value of the air valve setting should be close to the value that the PI controller finally settles down to. This vale can be read out under **Control PWM value (%)** and should then be entered to **offset%**.

In order to save these values to the non volatile memory of the board, press **Apply Changes**.

The spin can be tested from the web page using the start/stop toggle function.

VTU
Peripheral Control Board

Spin:

Spin frequenz: Start Stop

Set Spin frequenz[Hz] **target spin frequency**

Get Spin frequenz **actual spin frequency**

Get Spin Voltage **voltage of spin sensor**

Spin DC Voltage Counter	<input type="text" value="10"/>	safety thershold counter
Sample Down Threshold Voltage	<input type="text" value="0.35"/>	

Control: PWM Value [%] **current % air valve setting**

Offset[%]	<input type="text" value="30.99"/>	initial valve setting when spin starts	
P-factor	<input type="text" value="20.00"/>		proportional gain
I-factor	<input type="text" value="2.00"/>		integral gain
PostDelay in Ticks of time constant	<input type="text" value="3"/>		delay before control starts
time constant On Tick = 32ms	<input type="text" value="30"/>		

7.4.6.3 Setting the Spin Safety Threshold

The spin has a safety function that shuts the spin off as soon as the Lift is turned on either manually or during the initialization of the sample changer. In order to do this, the analog sensor that monitors the spinner has to be calibrated. It has three voltages for the three possible states:

No spinner	(typically 0.1 to 0.6V)
Black	(typically 0.6 to 2.5V)
White	(typically 2.5 to 4.5 V)

Since the values vary between shim stacks and also depends somewhat on the spinner type, it is necessary to calibrate this value. Usually this only has to be done once at installation.

1. Remove sample or activate the sample lift and note down the "empty" Spin Voltage.
2. Now insert or lower the sample and switch the spin on (either from the BSMS panel or from the web server)
3. Repeatedly read the Spin Voltage and note down the two values that can be read while the sample is spinning.
4. Now set the **sample down threshold voltage** to about half way between the lower of the two **black/white** spin voltages, and the **empty** spin voltage.

You can test if everything is set correctly. The spin should not cut off during normal operation, but should cut off immediately when the lift is activated.



The **Spin DC Voltage Counter** sets the number of times the value has to be below the threshold, before the safety shut off is activated. It should be set to 10 to avoid the spin cutting off on a single bad reading.



When the spin safety has been triggered, the spin has to be restarted manually, when the lift air has been turned off.

7.4.6.4 Calibrating and Measuring the Helium Level

When the measure button is pressed, the system lets a fixed measuring current flow through the measurement sensor for a time given by runtime x 32ms. The default is 100 x 32ms, i.e. 3.2 seconds.

After this time has elapsed, you can press **get value**. If the system has been properly calibrated, this will show the helium fill level in %

VTU

Peripheral Control Board

Helium Level Sensor:

runtime in ticks[32ms],
max Value = 200

Helium Level[%]

Helium Level Calibration:

- 1.) Insert the calibration resistor
- 2.) Press measure to read value
- 3.) Press corresponding get value Button

Helium: Full Level
[Voltage]

Helium: Low Level
[Voltage]

To calibrate the Helium Level Sensor :

1. Insert the **full** level calibration standard.
2. Click **measure**.
3. Wait minimum for the required time (see above).
4. Press **get value** at the corresponding field. This shows the voltage across the measurement resistor. For **full**, this should be **zero volts** (except for ECL 00 Peripheral Control Boards, in which case this value will be about **5**).
5. Repeat using the **empty** calibration standard.

6. For **empty**, this should be around **6 volts** (except for ECL 00 Peripheral Control Boards, in which case this value will be about **28**).
7. Press **Write to NVRAM**.

Helium Level Correction (ECL00 Peripheral Control Boards only)

ECL00 Peripheral Control Boards require an offset correction of 7%. This is entered in the section below the helium level calibration:

Helium Level Correction

factor offset

1.00 7.00 Set

For later boards*, this is not necessary and the setting should be:

Helium Level Correction

factor offset

1.00 0.00 Set

* in order to find out which version of board your system has, check the Firmware level of the Peripheral Control Board. V21 or lower is ECL00 and requires the calibration correction; all higher versions do not need the correction.

Function	Parameter	Method
BeAlive	Be alive	Get
Version	0.020	Get
LED	0x00	Set

7.4.6.5 Monitoring the Console and Shim Stack Temperature

VTU
Peripheral Control Board

Temperature:

Board temperature 31.00 Go


Shim temperature 22.02 Go

7.5 Additional Test Applications

The Fourier has some additional test applications the route certain RF signals back to the receiver in order to check if all the RF paths are working. These tests are based on the minispec software which is routinely installed with Fourier PCs.


The minispec Fourier Test Applications are found under:

C:\program files (x86)\Bruker the minispec\service_only\fourier test applications




FID.app

acquires an fid independently from topspin. Pulse length & frequency etc can be adjusted under **Settings**. TD, SWH are fixed for this test.




Fourier modulator test.app

Feeds the modulator signal back to the receiver and tests if a pulse is being generated.



The measurement is repeated 20 times and should give stable intensity results. This can be done both for 1H and for the C13 channel. Select the right channel by changing the **observe nucleus** in the window **Instrument Settings**.




Fourier modulator test.app

Feeds a small amount of the transmitter signal back to the receiver and tests if a pulse is being amplified. You can see the difference to the modulator signal by the rise time at the beginning. The signal intensity should be similar, but due to the variations of components cannot really be used to accurately measure the power.

The measurement is repeated 20 times and should give stable intensity results. This can be done both for 1H and for the C13 channel. Select the right channel by changing the **observe nucleus** in the window **Instrument Settings**.




Some oscillations on the top of the signal are normal and are caused by alias effects.



Fourier wobble.app

Since the **wobb** routine does not work with the Fourier, this app creates a similar wobble curve showing the resonance of the probe:




Configuration



The wobble routine is not very accurate and is NOT intended for adjusting the tune and match on the probe. This is factory adjusted and needs no adjustment in the field. The wobble routine is intended more for a quick check to make sure that cable, preamps and probe are functioning.

Fourier match.app

Gives an indication of the reflected power at a given frequency:



The values on the scale are arbitrary and do not represent the reflected power in %.

8 Pulse Programming

8.1 General Principle

The complete pulse programmer of the Fourier is implemented in the FPGA of the main board. This imposes some restriction on the variability and size of the pulse programs for this system. For example, there are no possibilities for branching in the pulse program. The pulse program can start due to an external trigger but once it started, it is not possible to halt/resume the program by to external triggers like the MAS signal. Another restriction is that RF phases can only be switched by 90° in a running experiment. This limits the phase cycles to the values, 0, 1, 2 and 3.

Unfortunately, it was not possible to use the output of the Topspin Avance compiler as an input for the Fourier pulse programmer due to a completely different topology. Therefore a different approach had to be used.

In addition to that, the console must also work with the miniSpec software because the same hardware is used as a high end miniSpec.

The solution is that the pulse program language for the Fourier is XML based. Instead of using the output of the Avance pulse program, only the pulse program parameters are extracted and an internal pulse program on the Fourier main board is started. As a consequence, for each Avance pulse program there must be a Fourier XML pulse program which generates exactly the same RF output. Both programs must have the same name. Starting the Avance pulse program loads the respective Fourier pulse program into the Fourier pulse programmer, sends the parameters like delays, pulse lengths and amplitudes, spectrum width, etc., to it and finally starts the experiment. The FIDs are sent to the workstation via Ethernet using the same CORBA protocol as the DRU on the Avance system.

The advantage is that the complete infrastructure of TopSpin can be used like used, NMRSIM, graphical pulse program display, the complete parameter handler, Cortab and PowerCheck. The operator works in the TopSpin environment only.

The disadvantage is that two pulse programs have to be written which must be perfectly synchronized. Any change on the Avance pulse program only takes effect on the experiment if the Fourier pulse program is changed as well. Bruker delivers a set of well tested pulse programs for the Fourier to cover the intended use of the Fourier spectrometer.

8.2 Structure of Pulse Programs

The Avance pulse program is written in a sequential form. Each delay or pulse on one channel influences the timing on the other channels. To illustrate this, the first pulses of an HMBC pulse program are used as an example.

Avance	Fourier	
All channels	Channel f1	Channel f12
d1	d1	d1
p1 ph1	p1 ph1	d31
d6	d6	d6
p3:f2 ph3	d31	p3 ph3
d0	d0	d0
p16:gp1	d36	p16:gp1
d16	d16	d16

Avance	Fourier	
p2 ph2	p2 ph2	d12
...
..

Table 8.1: Structure of Pulse Programs

with the definitions

d13=p1

d31=p3

d36=p16

To calculate the right timing for p3 on channel f2, e.g., the Avance pulse programmer must evaluate all pulses on channel f1 as well in addition to all delays. On the right-handed side, the parallel structure of the Fourier pulse program is shown. Both channels are completely independent of each other. The only synchronization is done at the beginning of the sequence. In order to achieve synchronicity throughout the complete sequence, the operator has to take care of the timing on both channels. E.g., a pulse on only one channel must be compensated by a delay of identical length in the other channel. As mentioned before, TopSpin still compiles the corresponding Avance style pulse program in order to generate all parameters correctly. As can be seen in the example above, additional delays are necessary for the synchronization of the Fourier pulse program channels. They appear as calculations in the TopSpin pulse programs.

In the pulse program text file, the complete code for the first channel is followed by that of the second one.

8.3 Syntax and Examples

The XML format used tags to assign code to a special function. There is always an opening and a closing tag. The code seems to be a bit complicated because all information for a pulse, gradient pulse or delay is completely defined.

Delays:

```
<DELAY>
  <Duration>
    <Unit>ms</Unit>
    <VarType>7</VarType>
    <JcampKwd>##$D[1]</JcampKwd>
  </Duration>
</DELAY>
```

As can be seen, a tag <DELAY> defines that the following information is a delay. The tag </DELAY> closes the definition. Usually, there are sub-definitions as the physical unit for the delay, e.g. Again, there is an opening tag, the information content and a closing tag. The tag <JcampKwd> indicates that this parameter is delivered by TopSpin after starting the experiment. In this case the parameter is d1.

Incremented delays:

```
<DELAY>
  <Duration>
```

```

    <Unit>ms</Unit>
    <VarType>7</VarType>
    <JcampKwd>##$D[0]+##$INF[1]</JcampKwd>
  </Duration>
</DELAY>

```

2D experiments typically increment an evolution time from slice to slice of the 2D. The syntax is that of a normal delay. The duration itself is calculated from slice to slice by adding the delay inf1.

Pulses:

```

<HF_PULSE>
  <Duration>
    <Unit>us</Unit>
    <VarType>7</VarType>
    <JcampKwd>##$P[1]</JcampKwd>
  </Duration>
  <Phase>
    <Unit>quadrant</Unit>
    <VarType>0</VarType>
    <Quantity>4</Quantity>
    <Val>0</Val>
    <Val>2</Val>
    <Val>1</Val>
    <Val>3</Val>
  </Phase>
  <Atten>
    <Unit>dB</Unit>
    <VarType>7</VarType>
    <JcampKwd>##$PL[1]</JcampKwd>
  </Atten>
  <Shape>
    <Unit>arb</Unit>
    <VarType>7</VarType>
    <Quantity>0</Quantity>
  </Shape>
</HF_PULSE>

```

The first tag defines the pulse duration. The next information is the phase cycle. The phase cycle is not defined at the end of the pulse program but is an integral part of the pulse definition. Please note, that the phases are written explicitly because they are fixed values for a given pulse program. So, they don't have to be extracted from the Jcamp file. The same syntax is also possible for durations or amplitudes by replacing the <JcampKw> tag by <VALUE> tags. The following information is the pulse amplitude. Once again, it is taken from

the Avance pre-compiler. In this case pl1. The last information is about the pulse shape. If it is a hard pulse, the shape quantity is "0", meaning that no information is available. Otherwise, the name of the shape file has to be given here.

If the phase cycle has to be incremented in a 2D experiment, an <Inc> tag has to be used. In the following example, the phase list is incremented by 1 for each new slice of the 2D experiment.

```
<Phase>
  <Unit>quadrant</Unit>
  <VarType>0</VarType>
  <Quantity>8</Quantity>
  <Val>0</Val>
  <Val>2</Val>
  <Val>1</Val>
  <Val>3</Val>
  <Inc>1</Inc>
</Phase>
```

Gradient pulses:

```
<GRAD_PULSE>
  <Direction>
    <Unit>number</Unit>
    <VarType>4</VarType>
    <Quantity>1</Quantity>
    <Val>4</Val>
  </Direction>
  <Duration>
    <Unit>us</Unit>
    <VarType>7</VarType>
    <JcampKwd>##$P[16]</JcampKwd>
  </Duration>
  <Gain>
    <Unit>%</Unit>
    <VarType>7</VarType>
    <JcampKwd>##$GPZ[1]</JcampKwd>
  </Gain>
  <Shape>
    <Unit>arb</Unit>
    <VarType>7</VarType>
    <Quantity>0</Quantity>
  </Shape>
</GRAD_PULSE>
```

The gradient pulse is defined by its **direction**, i.e., if it is a y, y or z gradient. It has a duration, p16 in this case, and an amplitude or gain of gpz1. It is a rectangular pulse, therefore, the shape name is 0.

Acquisition:

```

<RESET_ABS>
  <Execute> </Execute>
</RESET_ABS>

<ACQ_MULTI>
  <Duration>
    <Unit>ms</Unit>
    <VarType>7</VarType>
    <JcampKwd>##$AQ</JcampKwd>
  </Duration>
  <Phase>
    <Unit>quadrant</Unit>
    <VarType>0</VarType>
    <Quantity>4</Quantity>
    <Val>0</Val>
    <Val>2</Val>
    <Val>1</Val>
    <Val>3</Val>
  </Phase>
  <DataPoints>
    <Unit>points</Unit>
    <VarType>0</VarType>
    <JcampKwd>##$TD</JcampKwd>
  </DataPoints>
</ACQ_MULTI>

```

<RESETABS> switches to LO and starts the digitizer system. The timing of the acquisition is defined by <ACQ_MULTI>.

The acquisition duration is aq and it uses the phase cycle 0,2,1,3. TD data points are acquired. A dwell time or spectrum width is not necessary, they are derived from AQ and TD.

Loops:

```

<P_LOOP>
  <Start> </Start>
  <Counter>
    <Unit>counts</Unit>
    <VarType>4</VarType>
    <JcampKwd>##$L[4]</JcampKwd>

```

```
    </Counter>  
</P_LOOP>
```

....basic cycle....

```
<P_LOOP>  
    <Stop> </Stop>  
</P_LOOP>
```

Everything between the first and the second `<P_LOOP>` structure is repeated 14 times.

Channel assignment:

The complete 2 channel pulse program needs to know which pulses belong to which channel. Therefore, tags to define one complete channel are also necessary. The frequencies and transmitter paths for the channels are defined by edasp. The first channel is always f1, i.e., the acquisition channel.

```
<?xml version="1.0"?>
```

```
<PlsSetup>
```

```
<pulses>
```

```
code for channel f1....
```

```
</pulses>
```

```
<pulses>
```

```
code for channel f1....
```

```
</pulses>
```

```
</PlsSetup>
```

`<PlsSetup> </PlsSetup>` declares the text in between to be the pulse program.


`<pulses> </pulses>` defines the text in between to belong to one channel.

9 Wiring Schedules

Wiring Schedules


V1.4

Fourier 300 Verdrahtungsplan 1



V.1.4


Fourier_300_Verdrahtungsplan_2




Wiring Schedules

V1.4

Kourier 300 Verdrahtungsplan 3



V.1.4



FOURIER 300
KONSOLE EINHEIT 300MHZ
(E1451100)
KLEMMENPLAN

21.02.2011

LINKE SEITE		FRONTANSICHT		RECHTE SEITE	
SAP-KOMPONENTENBEZEICHNUNG	SAP-TEILNUMMER	DC-VERTEILER- KLEMMLEISTE	SAP-KABELBEZEICHNUNG	SAP-KABELBEZEICHNUNG	SAP-KOMPONENTENBEZEICHNUNG
SAP-TEILNUMMER	SAP-TEILNUMMER	SAP-TEILNUMMER	SAP-TEILNUMMER	SAP-TEILNUMMER	SAP-TEILNUMMER
RESERVE	RESERVE	PE ₁	RESERVE	RESERVE	RESERVE
RESERVE	RESERVE	PE ₂	RESERVE	RESERVE	RESERVE
RESERVE	RESERVE	PE ₃	RESERVE	RESERVE	RESERVE
SW MODEP. S. 500W 22-30V 18.5A 89004	(g/n/9)	PE ₄	PE 1*0,75mm ²	DSPEC PNEUMATIKVERTEILERBLOCK + DIGITAL 1804214	
VENT SCHUTZSTIFTER LZ 23 (11094)	KABEL SPEZ GND 300 4-ADRIG E1451233	PE ₅	PE 1*0,75mm ²	DSPEC PNEUMATIKVERTEILERBLOCK + DIGITAL 1804214	
DSPEC KONSOLE RUECKWAND (E1451165)	KABEL SPEZ GND 300 4-ADRIG E1451233	PE ₆	KABEL SPEZ GND 8607-ADRIG E1451234	DSPEC KONSOLE BOOENPLATTE (E1451166)	
ST JEC APPLIANCE PLUG BU SHUTTERED (67902)	KABEL SPEZ GND 860 6-ADRIG E1451235	PE ₇	KABEL SPEZ GND 8607-ADRIG E1451234	DSPEC VENTILBLOCK HALTEWAHSEL (E1451157)	
SHIM STEUERBOARD TRAEGERPLATTE (E1451125)	KABEL SPEZ GND 860 6-ADRIG E1451235	PE ₈	KABEL SPEZ GND 8607-ADRIG E1451234	KONSOLE MICRO HUB BEFESTIGUNG (E7002124)	
LED STEUER-80 BEFESTIGUNG (E7001096)	KABEL SPEZ GND 860 6-ADRIG E1451235	PE ₉	KABEL SPEZ GND 8607-ADRIG E1451234	SIGNAL CHAINWEL BOOENPLATTE (E7003313)	
SENDERBOARD TRAEGERPLATTE (E1451127)	KABEL SPEZ GND 860 6-ADRIG E1451235		KABEL SPEZ GND 8607-ADRIG E1451234	STEUERBOARD TRAEGERPLATTE (E1451128)	
LOCKBOARD TRAEGERPLATTE (E1451128)	KABEL SPEZ GND 860 6-ADRIG E1451235			DSPEC KONSOLE STUEZSTREBE 3 (E1451168)	
DSPEC KONSOLE STUEZSTREBE (E1451161)				DSPEC KONSOLE STUEZSTREBE 4 (E1451164)	
DSPEC KONSOLE STUEZSTREBE 2 (E1451162)					

FOURER 300
KONSOLE EINHEIT 300MHZ
(E1451100)
KLEINHEBPLAN

21.02.2011

LINKE SEITE		FRONTANSICHT		RECHTE SEITE	
SAP-KOMPONENTENBEZEICHNUNG	SAP-KABELBEZEICHNUNG	DC-VERTEILER- KLEINLEISTE	SAP-KABELBEZEICHNUNG	SAP-KOMPONENTENBEZEICHNUNG	SAP-TEILNUMMER
SCHALTNETZTEIL -45+12-12+24V 5/0/0/0A	KABEL SPEZ 600 U-VERTEILER > P.S.	•	KABEL SPEZ 600 U-VERT > VTU	VT STEUERBOARD	E1451222
SCHALTNETZTEIL -45+12-12+24V 5/0/0/0A	KABEL SPEZ 600 U-VERTEILER > P.S.	•	KABEL SPEZ 600 U-VERT > VTU	VT STEUERBOARD	E1451222
SCHALTNETZTEIL -45+12-12+24V 5/0/0/0A	KABEL SPEZ 600 U-VERTEILER > P.S.	•	KABEL SPEZ 600 U-VERT > VTU	VT STEUERBOARD	E1451222
VENT 24VDC 6AC/5MH 92X92/25	schwarz GND 1*0,25mm ²	•	schwarz GND 1*0,25mm ²	VENT 24VDC 6AC/5MH 92X92/25	E1451225
DSPEC X.VORVERSTAEKER (13C)	KABEL SPEZ 600 +12V UVERT > 13C VORVERST.	•	(schwarz/rot/blau/rot) GND 1 * 0,25mm ²	ETHERNET SWITCH 10/100 MBIT 8RJ45	E1451226
DSPEC VORVERSTAEKER 300MHZ	KABEL SPEZ 750 +12V UVERT > H VORVERST.	•	KABEL SPEZ 550 +5V UVERT > GA8 BACKPL	DSPEC GA82 BACKPLANE	E1451226
DSPEC TASTER STEUERBOARD	KABEL SPEZ 600 LED CTRL > P.S. SV	•	KABEL SPEZ 550 +5V UVERT > GA8 BACKPL	DSPEC GA82 BACKPLANE	E1451226
DSPEC PERIPHERIE-STEUERBOARD	KABEL SPEZ 600 2428V UVERT > PERIPH. BD	•	KABEL SPEZ 900 +12V UVERT > GA8 BACKPL	DSPEC GA82 BACKPLANE	E1451226
DSPEC LOCKSCHALTERBOARD	KABEL SPEZ 710 DUBOX BU8U UVERT > BOX	•	RESERVE	RESERVE	E1451225
RESERVE	RESERVE	•	RESERVE	RESERVE	E1451225
SCHALTNETZTEIL -45+12-12+24V 5/0/0/0A	KABEL SPEZ 600 U-VERTEILER > P.S.	•	KABEL SPEZ 900 U-VERT > VTU	VT STEUERBOARD	E1451222
DSPEC TASTER STEUERBOARD	KABEL SPEZ 600 LED CTRL > P.S. SV	•	KABEL SPEZ 550 +5V UVERT > GA8 BACKPL	DSPEC GA82 BACKPLANE	E1451226
RESERVE	RESERVE	•	KABEL SPEZ 550 +5V UVERT > GA8 BACKPL	DSPEC GA82 BACKPLANE	E1451226
SCHALTNETZTEIL -45+12-12+24V 5/0/0/0A	KABEL SPEZ 600 U-VERTEILER > P.S.	•	KABEL SPEZ 900 U-VERT > VTU	VT STEUERBOARD	E1451222
DSPEC X.VORVERSTAEKER (13C)	KABEL SPEZ 600 +12V UVERT > 13C VORVERST.	•	KABEL SPEZ 550 +5V UVERT > GA8 BACKPL	DSPEC GA82 BACKPLANE	E1451226
DSPEC VORVERSTAEKER 300MHZ	KABEL SPEZ 750 +12V UVERT > H VORVERST.	•	(schwarz/rot/blau/rot) +12V 1 * 0,25mm ²	ETHERNET SWITCH 10/100 MBIT 8RJ45	E1451225
DSPEC LOCKSCHALTERBOARD	KABEL SPEZ 710 DUBOX BU8U UVERT > BOX	•	RESERVE	RESERVE	E1451225
SCHALTNETZTEIL -45+12-12+24V 5/0/0/0A	KABEL SPEZ 600 U-VERTEILER > P.S.	•	RESERVE	RESERVE	E1451225
VENT 24VDC 6AC/5MH 92X92/25	parallel 3*2400 zwischen Kabel und Klemmleiste rot +24V 1*0,25mm ²	•	parallel 3*2400 zwischen Kabel und Klemmleiste rot +24V 1*0,25mm ²	VENT 24VDC 6AC/5MH 92X92/25	E1451225
DSPEC PERIPHERIE-STEUERBOARD	KABEL SPEZ 600 2428V UVERT > PERIPH. BD	•	KABEL SPEZ 900 U-VERT > VTU	VT STEUERBOARD	E1451222
SCHALTNETZTEIL -45+12-12+24V 5/0/0/0A	KABEL SPEZ 600 U-VERTEILER > P.S.	•	RESERVE	RESERVE	E1451222
RESERVE	RESERVE	•	RESERVE	RESERVE	E1451222
RESERVE	RESERVE	•	RESERVE	RESERVE	E1451222

2 von 3

FOURIER 300
 KONSOLE EINHEIT 300MHZ
 (E1451100)
 KLEMMENPLAN
 21.02.2011

LINKE SEITE		FRONTANSICHT		RECHTE SEITE	
SAP-KOMPONENTENBEZEICHNUNG	SAP-KABELBEZEICHNUNG	DC-VERTEILER- KLEMMENLISTE	SAP-KABELBEZEICHNUNG	SAP-KOMPONENTENBEZEICHNUNG	SAP-TEILENUMMER
SAP-TEILENUMMER	SAP-TEILENUMMER		SAP-TEILENUMMER	SAP-TEILENUMMER	
SW MODE P. S. 500W 22-30V 18.5A 89064	KABEL SPEZ 910 *28V U-VERTEILER > SMC E1451213	+ 3M 24 ●	KABEL SPEZ 900 *28V U-VERT > VTU BD E1451223	VT STEUERBOARD E7003520	
SW MODE P. S. 500W 22-30V 18.5A 89064	KABEL SPEZ 910 *28V U-VERTEILER > SMC E1451213	+ 3M 25 ●	KABEL SPEZ 900 *28V U-VERT > VTU BD E1451223	VT STEUERBOARD E7003520	
SW MODE P. S. 500W 22-30V 18.5A 89064	KABEL SPEZ 910 *28V U-VERTEILER > SMC E1451213	+ 3M 26 ●	KABEL SPEZ 730 *28V U-VERT > DSPEC BD E1451214	DSPEC STEUERBOARD E1451110	
SHM SPULE 5V-STEUER-BOARD E1423120	KABEL SPEZ 330 *28V U-VERTEIL > SHM BD E1451216	+ 3M 28 ●	KABEL SPEZ 730 *28V U-VERT > DSPEC BD E1451214	DSPEC STEUERBOARD E1451110	
SHM SPULE 5V-STEUER-BOARD E1423120	KABEL SPEZ 330 *28V U-VERTEIL > SHM BD E1451216	+ 3M 29 ●	KABEL SPEZ 540 *28V U-VERT > GA.B BACKPL. E1451224	DSPEC GA.B2 BACKPLANE E1451170	
DSPEC LOCK CONTROL BOARD E1451180	KABEL SPEZ 530 *28V U-VERTEIL > LOCK BD E1451215	+ 3M 30 ●	RESERVE	RESERVE	
DSPEC LOCK CONTROL BOARD E1451180	KABEL SPEZ 530 *28V U-VERTEIL > LOCK BD E1451215	+ 3M 31 ●	RESERVE	RESERVE	
DSPEC PERIPHERIE-STEUERBOARD E1451130	KABEL SPEZ 800 *28V U-VERT > PERIPH. BD E1451217	+ 3M 32 ●	RESERVE	RESERVE	
RESERVE	RESERVE	+ 3M 33 ●	RESERVE	RESERVE	
SW MODE P. S. 500W 22-30V 18.5A 89064	KABEL SPEZ 910 *28V U-VERTEILER > SMC E1451213	OND 34 ●	KABEL SPEZ 900 *28V U-VERT > VTU BD E1451223	VT STEUERBOARD E7003520	
SW MODE P. S. 500W 22-30V 18.5A 89064	KABEL SPEZ 910 *28V U-VERTEILER > SMC E1451213	OND 35 ●	KABEL SPEZ 900 *28V U-VERT > VTU BD E1451223	VT STEUERBOARD E7003520	
SW MODE P. S. 500W 22-30V 18.5A 89064	KABEL SPEZ 910 *28V U-VERTEILER > SMC E1451213	OND 36 ●	KABEL SPEZ 730 *28V U-VERT > DSPEC BD E1451214	DSPEC STEUERBOARD E1451110	
SW MODE P. S. 500W 22-30V 18.5A 89064	KABEL SPEZ 910 *28V U-VERTEILER > SMC E1451213	OND 37 ●	KABEL SPEZ 730 *28V U-VERT > DSPEC BD E1451214	DSPEC STEUERBOARD E1451110	
SHM SPULE 5V-STEUER-BOARD E1423120	KABEL SPEZ 330 *28V U-VERTEIL > SHM BD E1451216	OND 38 ●	KABEL SPEZ 540 *28V U-VERT > GA.B BACKPL. E1451224	DSPEC GA.B2 BACKPLANE E1451170	
SHM SPULE 5V-STEUER-BOARD E1423120	KABEL SPEZ 330 *28V U-VERTEIL > SHM BD E1451216	OND 39 ●	RESERVE	RESERVE	
DSPEC LOCK CONTROL BOARD E1451180	KABEL SPEZ 530 *28V U-VERTEIL > LOCK BD E1451215	OND 40 ●	RESERVE	RESERVE	
DSPEC LOCK CONTROL BOARD E1451180	KABEL SPEZ 530 *28V U-VERTEIL > LOCK BD E1451215	OND 41 ●	RESERVE	RESERVE	
DSPEC PERIPHERIE-STEUERBOARD E1451130	KABEL SPEZ 800 *28V U-VERT > PERIPH. BD E1451217	OND 42 ●	RESERVE	RESERVE	
RESERVE	RESERVE	OND 43 ●	RESERVE	RESERVE	

10 Maintenance and Repair

The Console is almost maintenance free. It is recommended to regularly (every 6 to 12 months) depending on the environment to check and if necessary clean the filter mat for the ventilators. The filter tray, see the next figure (1), can be removed without taking the console cover off.




Figure 10.1: The filter tray

The filter itself can be cleaned with ordinary tap water and should be dried before replacing. Do not operate the instrument without filter for prolonged times, since dirt will build up in the electronics. Repairs on the console should only be carried out by trained Bruker personnel. For probe and magnet maintenance, refer to the relevant manuals.

11 Cleaning and Disposal

11.1 Cleaning

Do not use any detergent or other cleaning solvents. Use only water or neutral cleaning fluids. Usage of volatile cleaners like thinner or benzine may damage the surface of the unit.

- Clean the outside of the device with a soft, lint-free cloth dampened in water.



Wait until the unit is completely dry before you reconnect the power cable.

11.2 Disposal

Following the end of its useful life, the device must be dismantled and disposed of in accordance with the environmental regulations.



Installation, initial commissioning, retrofitting, repairs, adjustments or dismantling of the device must only be carried out by employees of the manufacturer or persons authorized by the manufacturer.

11.3 Safety

Electrical System

WARNING



Danger of injury from electrical shock!

A life threatening shock may result when the housing is open during operation.

- ▶ Only qualified personnel should open the housing.
- ▶ Disconnect the device from the electrical power supply before opening the device. Use a voltmeter to verify that the device is not under power!
- ▶ Be sure that the power supply cannot be reconnected without notice.

Improper Dismantling



WARNING

Danger of injury due to improper dismantling!

Stored residual energy, angular components, points and edges on and in the device or on the tools needed can cause injuries.

- ▶ Ensure sufficient space before starting work.
- ▶ Handle exposed, sharp-edged components with care.
- ▶ Dismantle the components properly.
- ▶ Secure components so that they cannot fall down or topple over.
- ▶ Consult the manufacturer if in doubt.

11.4 Dismantling

Before starting dismantling:

- Shut down the device and secure to prevent restarting.
- Physically disconnect the power supply from the device; discharge stored residual energy.
- Remove consumables, auxiliary materials and other processing materials and dispose of in accordance with the environmental regulations.
- Dismantle the device by following the installation instructions in reverse.

Clean assemblies and parts properly and dismantle in compliance with applicable local occupational safety and environmental protection regulations.

11.5 Disposal Instructions

If no return or disposal agreement has been made, send the dismantled components for recycling.

- Scrap metals.
- Send plastic elements for recycling.
- Sort and dispose of other components in accordance with their material composition.

NOTICE

Danger to the environment from incorrect handling of pollutants!

Incorrect handling of pollutants, particularly incorrect waste disposal, may cause serious damage to the environment.

- ▶ Always observe local environmental regulations regarding handling and disposal of pollutants.
- ▶ Take the appropriate actions immediately if pollutants escape accidentally into the environment. If in doubt, inform the responsible municipal authorities about the damage and ask about the appropriate actions to be taken.

12 Technical Specifications

Data	Value	Unit
Voltage	110-230	V AC
Apparent power consumption, maximum	300	VA
Circuit protection	2 x 4AT 250V	A
Frequency	50/60	Hz

Table 12.1: Electrical Connection Values

Transmit Channels	
Frequency Range	5 – 405 MHz
No of Channels	2
Frequency Resolution	0.25 Hz
No of Phases	4
Timing Resolution	20 ns
Transmit Power for 1H	ca. 50 W
Transmit Power for X Channel (13C, 15N)	ca. 200W
Power Level Adjustment	0 - 60 dB in 1dB steps
Amplitude Variation	14-bit

Table 12.2: Transmit Channels

Receiver Channel	
No of Channels	1
ADC	50 MHz 14-bit
Effective Resolution	16.5 bits per complex pair at 1 MHz
Detection Type	Direct Digital Quad Detection
No of Phases	4
Analog Bandwidth	3 MHz
Digital Filter	FIR Filter 1kHz to 6 MHz in steps of factor 2
Max No of points per acquisition	2 x 32k (after decimation)

Table 12.3: Receiver Channel

Gradient Control	
No. of Channels	1 internal (GAB/2), 3 external possible
Gradient Interface	LVDS (80 MHz)

Table 12.4: Gradient Control

Shim System	
No of Shim Channels	20 (BOSS1 compatible)
Shim Resolution	16-bit
Shim Current	max 300 mA per channel
Shim Temperature Monitor	yes

Table 12.5: Shim System

Lock System	
Lock Nucleus	2H
Lock Type	digital P+I
Lock Speed	1 kHz
Lock Power	-30dBm to 0 dBm in 1 dB steps
Lock Resolution	0.025 Hz
Lock Range	10 kHz

Table 12.6: Lock System

Peripheral Control	
Helium Level Sensor	yes
Sample Spin Control	yes
Sample Lift Control	yes
ATMA Control	no

Table 12.7: Peripheral Control

N2 Temperature Controller	
Temperature Control	Digital PID
Temperature Range	100 to 500 °K
Temperature Sensor	PT100 or Type K Thermocouple
Temperature Stability	0.1 °
Gas Flow Regulation	0 to 1600 l/h
N2 Evaporator Control	yes

Table 12.8: N2 Temperature Controller

Software & Environment	
Spectrometer Control Software	TopSpin, minispec software
Interface	100 MB/s Ethernet
Embedded Web Server	yes

Software & Environment	
Pulse Programmes	predefined, with TopSpin variable parameters
1D Experiments	ZG, ZG30, ZGIG(1,2), ZGPG(2), ZGPG30(2) DECP90, APT, APTJR, DEPT (1) GARP & GARP4 decoupling, (2) WALTZ64 decoupling
2D Experiments	COSYQF, COSYGPQF, COSYPH TOCSYPH, NOESYPH, NOESYGPPH HSQCGPPH(1), HMQCGPQF(1), HMBCGPNDQF(1), HMBCGPLNDQF(1)

Table 12.9: Software and Environment

13 Contact

Manufacturer:

Bruker BioSpin NMR
Silberstreifen
D-76287 Rheinstetten
Germany
Phone: +49 721-5161-6155
<http://www.bruker.com>
WEEE DE43181702

NMR Hotlines

Contact our NMR service centers.

Bruker BioSpin NMR provide dedicated hotlines and service centers, so that our specialists can respond as quickly as possible to all your service requests, applications questions, software or technical needs.

Please select the NMR service center or hotline you wish to contact from our list available at:

<http://www.bruker.com/service/information-communication/helpdesk/magnetic-resonance.html>

List of Figures

Figure 2.1:	The Fourier 300 System	11
Figure 5.1:	The Fourier 300 Console	27
Figure 5.2:	Fourier 300 Console Overview	28
Figure 5.3:	Console: Left Side	29
Figure 5.4:	Console: Right Side	30
Figure 5.5:	Console: Front View	31
Figure 5.6:	Location of the Console Air Filter Tray	32
Figure 5.7:	Console Input and Output Connections	33
Figure 5.8:	DDS2	34
Figure 5.9:	DDS1	35
Figure 5.10:	Two Shim Boards	36
Figure 5.11:	Lock Switch.....	39
Figure 5.12:	Lock Switch Block Diagram	40
Figure 5.13:	Lock Switch connections & LEDS.....	40
Figure 5.14:	Ventilator.....	41
Figure 5.15:	Connections.....	42
Figure 6.1:	Checking the DHCP Server	44
Figure 6.2:	cf.....	46
Figure 6.3:	POWCHK.....	47
Figure 6.4:	expinstall.....	50
Figure 6.5:	BSMS Panel	58
Figure 6.6:	Result of Flock with one Solvent Peak	61
Figure 6.7:	Gradient echo	65
Figure 7.1:	Fixed User Specific IP Address	69
Figure 7.2:	Power Up Jumper	70
Figure 7.3:	Wide Coarse too high	91
Figure 7.4:	Wide DAC too low.....	91
Figure 7.5:	Wide DAC OK.....	91
Figure 7.6:	Phase wrong.....	92
Figure 7.7:	Phase Correct.....	92
Figure 10.1:	The filter tray.....	119

List of Tables

Table 1.1: Font and Format Conventions	9
Table 2.1: Overview Installation and Operation Requirements for Personnel	12
Table 5.1: Preamplifiers	38
Table 6.1: Hard Specs for the Fourier	67
Table 6.2: Typical Values for the Fourier	68
Table 7.1: Status and Description	89
Table 8.1: Structure of Pulse Programs	105
Table 12.1: Electrical Connection Values	123
Table 12.2: Transmit Channels	123
Table 12.3: Receiver Channel	123
Table 12.4: Gradient Control	124
Table 12.5: Shim System	124
Table 12.6: Lock System	124
Table 12.7: Peripheral Control	124
Table 12.8: N2 Temperature Controller	124
Table 12.9: Software and Environment	124

Index

A

Accessories 42

B

BOSS1 Shim System 11

C

Can-Bus 36

cf 45

D

Dismantling 122

F

Fourier 300 11

Fourier 300 Console 27

Fourier 300 spectrometer 11

Fourier match.app 104

Fourier modulator test.app 102

Fourier wobble.app 103

G

GRADSHIM 47

L

Laboratory personnel 15

Lock Switch 38

Lock System 35

LVDS 37

M

Main Board 34

N

Nuclei List 47

P

Personal protective equipment 16

Personnel qualifications 15

Power Supplies 40

R

RF routing 49

S

Shim Control 36

Solvents 52

Spin DC Voltage Counter 97

T

tty Settings 46


V

Variable Temperature Unit 36

Ventilators 41

VTU 36





Bruker Corporation

info@bruker.com
www.bruker.com

Order No: Z31980